




PRODUCT CATALOGUE

PARAMOUNT
QUEST FOR BETTER HEALTHCARE



The Paramount Philosophy:

With a vision of being a global leader in the medical disposable industry, Paramount Surgimed believes in providing its customers a timely supply of quality products at a competitive price. The company came into existence in 1993 under the Grover Group of companies capitalizing on the rich and diverse experience of over 60 years. Now, the brand is symbolized with trust and quality in the medical industry. The company started with the manufacturing of Surgical Blades and is now a leading supplier of a huge range of disposable products.

Paramount exports to over 97 countries, has 40 distributors and 400 Dealers across India. Our State-of-the-Art production facility compares to any international standard and is accredited with ISO: 13485:2016, CE Certificate and US FDA Device Listing.

The Paramount Vision:

Paramount Surgimed and all its employees share the vision of their dynamic Managing Director, Mr. Shaily Grover to be a global leader in the Medical Industry. Mr. Grover has ensured that the Group's reputation amongst its associates is synonymous with top quality, updated technology, dedicated approach and customer friendliness.



Corporate Office – Okhla, New Delhi



Factory – Bhiwadi, Rajasthan



Factory – Bawal, Haryana



Factory – Bhiwadi, Rajasthan

Surgical Blades	02
Stitch Cutters	04
Thumb Scalpel	05
Dermal Blade	06
Gouge Blades & Handles	07
Disposable Scalpels	08
Safety Scalpel	09
Universal Safety Scalpel	10
Biopsy Punch	11
Dermal Curette	13
Ophthalmic Knives	15
Chisel/Fine/Micro Surgical Blade	17
Skin Graft Blade	18
Blood Lancet	18
Myringotomy Knife	19
Face Mask 3 PLY/ Type IIR	20
Face Mask N95/ FFP2	21
Disposable Adult Diapers	22
Under Pad - AroKleen & PupPee	23
Patient Wipes, Baby Wipes	24
Disposable Prep Razor	25
Lung Exerciser, Pulse Oximeter, Infrared Thermometer	26
Bed Bath Wipes, Baby Wipes	27
Cleansing Wipes, Germ Protection Wipes	28
Pet Wipes, Hot Water Bottle	29
Certificates	30

Surgical Blades

Carbon Steel & Stainless Steel

Features

- Stainless Steel Surgical Blades and Carbon Steel Surgical Blades.
- Packed in peel open pouches with a VCI protective inner paper.
- Shelf life of 5 years - Sterile Product
- Stainless Steel (Grade F as per ISO 7153-1:2016)
- Carbon Steel (Grade- As per BS 2982:1992)
- Complies the Hardness requirement of BS 2982:1992 standard.
- Complies the fitment requirement of ISO 7740:1985 standard requirement.



FITMENT 3



FITMENT 4





Details:

Product Name:	Surgical Blades
Size of Inner Box (cm):	8.5 x 7.0 x 2.7
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	40 x 39 x 16.5
Per 20 feet container:	825 Carton

Size	Carbon Steel	Stainless Steel	Fitment No
9	PCS9	PSS9	3
10	PCS10	PSS10	3
10A	PCS10A	PSS10A	3
10R	PCS10R	PSS10R	3
11	PCS11	PSS11	3
11K	PCS11K	PSS11K	3
11P	PCS11P	PSS11P	3
12	PCS12	PSS12	3
12D	PCS12D	PSS12D	3
13	PCS13	PSS13	3
14	PCS 14	PSS14	3
15	PCS15	PSS15	3
15B	PCS15B	PSS15B	3
15C	PCS15C	PSS15C	3
15D	PCS15D	PSS15D	3
15S	PCS15S	PSS15S	3
16	PCS16	PSS16	3
17	PCS17	PSS17	3
18	PCS18	PSS18	4
19	PCS19	PSS19	4
20	PCS20	PSS20	4
21	PCS21	PSS21	4
22	PCS22	PSS22	4
22A	PCS22A	PSS22A	4
23	PCS23	PSS23	4
24	PCS24	PSS24	4
24D	PCS24D	PSS24D	4
25	PCS25	PSS25	4
26	PCS26	PSS26	4
36	PCS36	PSS36	4
36D	PCS36D	PSS36D	4
40	PCS40	PSS40	-
40B	PCS40B	PSS40B	-
60	PCS60	PSS60	8
60B	PCS60B	PSS60B	8
70	PCS70	PSS70	8

Post Mortem Blades –Size 60

Features:

- Used in the fields of Pathology and Histopathology for autopsy
- It is used with handle No.8
- Available in both material carbon steel and stainless steel
- Box of 50
- Individual peel pouch



Post Mortem Blades –Size 40

Features:

- The PM40 range of blades are meant to be used with screw type heavy duty handles
- Heavy duty blade
- Box of 10



Stitch Cutters

Stainless Steel & Carbon Steel

Description

Stitch Cutter is used for general suture removal. Appearing similar to a scalpel, this instrument requires no handle and is a simple and convenient, way to remove basic sutures. The slightly curved, flat, bull nose tip slides under the suture and when rotated, gently removes the suture from the site. Long stitch cutter is made of stainless steel and short stitch cutter is made of carbon steel. The mini stitch cutter is made of stainless steel can be used with handle No.3L, 7 & 9. Alternatively, it is also available on pre-mounted ABS disposable plastic handle.

Features

- 5 Year Shelf Life – Sterile Product
- Each blade is protected by VCI paper in a peel-open foil sachet.
- Convenience of single-use disposables.
- Sizes clearly marked on each blade.
- Each box contains 100 pcs

Size	Carbon Steel	Stainless Steel	Handle No.
LSC	–	PLSCSS	-
SSC	PSSCCS	–	-
MSC	–	PMSCSS	PCS3L

Details:

Product Name:	Mini Stitch Cutter
Size of Inner Box (cm):	8.5 x 7.0 x 2.7
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	40 x 39 x 16.5
Per 20 feet container:	825 Carton



Details:

Product Name:	Short Stitch Cutter
Size of Inner Box (cm.):	11 x 8 x 3.5
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm.):	46.5 x 40.0 x 25.5
Per 20 feet container:	550 Carton



Details:

Product Name:	Long Stitch Cutter
Size of Inner Box (cm):	14.8 x 7 x 3.1
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	45 x 34 x 36
Per 20 feet container:	475 Carton

**Details:**

Product Name:	Mini Stitch Cutter with Handle
Size of Inner Box (cm):	14.5 x 70 x 3.5
Unit/Box:	Bulk Pack of 50 Pcs
Outer Boxes:	Not Using
Export carton:	10 Bags of 50 Pcs x 10 Outer
Export Carton size (cm):	61 x 50 x 3.5
Per 20 feet container:	290 Carton




Thumb Scalpel

Stainless Steel

Description

Disposable thumb scalpel or pocket scalpel is used in place of a surgical blade. This reduces or even prevents the incidences of injuries that can occur during the process of fitting or removing contaminated blades from handles.

The disposable thumb scalpels have a protective cap/guard for the blade that can be used so the scalpel can be safely disposed in the operation theater/room disposal container.



Dermal Blade


Stainless Steel blade

Description

- The device is intended for cutaneous surgery giving smooth clean excisions with less tissue trauma. It is Used to remove surface protuberances, flat lesions, or deep lesions by saucer-shaped incision. Dermal blade is made of stainless-steel blade mounted on flexible plastic holder. The device is under category of surgically invasive device and for transient use. It is covering under surface device (nature of body contact) and contact with intact skin and mucosal membrane. The device is invasive device and for transient use.

Details:


Product Name:	Dermal Biopsy Blade
Size of Inner Box (cm):	15.5 x 11.5 x 11
Unit/Box:	20 Pcs
Export carton:	20 Boxes
Export Carton size (cm):	59 x 32 x 22.5
Per 20 feet container:	700 Carton



Metal Handle

Description

- Metal handles comes in different sizes and are used according to the applications. The most popular stainless steel handle sizes used are size #3 and size #4 handles.



Gouge Blades & Handles

Carbon Steel & Stainless Steel

Description

These surgical blades are used for foot surgeries and podiatry purposes. Gouge Blades are used by various podiatrists for hand, feet and skin surgery. They are available in the following sizes: Nos. 1, 2, 3, 4, 5, 6, 8, 10, 12, 15 and 18. Each blade is protected by VCI paper in a peel-open foil sachet. Gouge blades fit into handle no. 3 and special gouge handles.



Details:

Product Name: Gouge Blades No. 1 to 5
Size of Inner Box (cm): 10 x 10 x 3.3
Unit/Box: 100 Pcs
Outer Boxes: 10 Boxes x 100 Units
Export carton: 10 Boxes x 10 Outer
Export Carton size (cm): 59.5 x 47 x 21
Per 20 feet container: 425 Carton

Size	Product Code	Product Code	Size	Product Code	Product Code	Fittment No	Handle
1	PGCS1	PGSS1	1V	PGCS1V	PGSS1V	3	PGMH
2	PGCS2	PGSS2	2V	PGCS2V	PGSS2V	3	PGMH
3	PGCS3	PGSS3	3V	PGCS3V	PGSS3V	3	PGMH
4	PGCS4	PGSS4	1R	PGCS1R	PGSS1R	3	PGMH
5	PGCS5	PGSS5	2R	PGCS2R	PGSS2R	3	PGMH
6	PGCS6	PGSS6	3R	PGCS3R	PGSS3R	3	PGMH
8	PGCS8	PGSS8				3	PGMH
10	PGCS10	PGSS10				3	PGMH
12	PGCS12	PGSS12				3	PGMH
15	PGCS15	PGSS15				3	PGMH
18	PGCS18	PGSS18				3	PGMH

Details:

Product Name: Gouge Blades No. 6 to 18
Size of Inner Box (cm): 13.8 x 10 x 3.3
Unit/Box: 100 Pcs
Outer Boxes: 10 Boxes x 100 Units
Export carton: 10 Boxes x 10 Outer
Export Carton size (cm): 57 x 39 x 30
Per 20 feet container: 390 Carton




Disposable Scalpels

Stainless Steel blade

Features

- Pre-assembled stainless steel surgical blades mounted on ABS plastic handle.
- The device is used to prevent injury during fixing the blade on the metal handle.
- Fluted finger grips and the blades are safely protected by LDPE guard.
- Individually packaged in a peel pouch.
- Convenience of single-use.
- Size clearly marked on each scalpel blade on blister pouch.
- Each dispenser box contains 10 scalpels.
- Shelf life of 5 years - Sterile Product



Details:

Product Name:	Disposable Scalpel
Size of Inner Box (cm):	21 x 7 x 3
Unit/Box:	10 Pcs
Export carton:	100 Boxes
Export Carton size (cm):	44 x 37 x 33
Per 20 feet container:	370 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
9	PDSS9	3	18	PDSS18	4
10	PDSS10	3	19	PDSS19	4
10A	PDSS10A	3	20	PDSS20	4
11	PDSS11	3	21	PDSS21	4
11K	PDSS11K	3	22	PDSS22	4
12	PDSS12	3	22A	PDSS22A	4
13	PDSS13	3	23	PDSS23	4
14	PDSS14	3	24	PDSS24	4
15	PDSS15	3	25	PDSS25	4
15B	PDSS15B	3	36	PDSS36	4
15C	PDSS15C	3			
16	PDSS16	3			
17	PDSS17	3			



Safety Scalpel

Stainless Steel blade


Features

- It is a safe alternative to the traditional scalpel and fulfills the requirement to reduce sharp injuries.
- Ribbed on both side of the sheath which provide extra gripping for left handed as well as right handed surgeon.
- One-handed activation allowing the surgeon to focus on the surgical site.
- Audible clicks makes the scalpel more secure and protective.
- "Used & Un-Used" indications prevent contamination and transmission of blood borne disease.
- Permanent safety lock prevents accidental cuts.
- Latex Free.
- Measurement scale 0 to 5 cm is embossed on the handle.
- Ergonomically designed safety scalpel.
- ABS plastic handle and sheath gives extra strength to the product.
- Universal size handle that is ideal for kits, trays & packs.

Details:

Product Name:	Safety Scalpel
Size of Inner Box (cm):	21 x 7 x 3
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	44 x 37 x 17
Per 20 feet container:	300 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
10	PNSF10	3	21	PNSF21	4
10A	PNSF10A	3	22	PNSF22	4
11	PNSF11	3	23	PNSF 23	4
11K	PNSF11K	3	24	PNSF24	4
13	PNSF13	3	25	PNSF25	4
15	PNSF15	3	36	PNSF36	4
15B	PNSF15B	3	MSC	PNSFMSC	3
15C	PNSF15C	3			
17	PNSF17	4			
18	PNSF18	4			
19	PNSF19	4			
20	PNSF20	4			



Universal Safety Scalpel

Stainless Steel blade

Features

- **Safety Mode:** When you retract back after fully extending the blade, it makes a click sound to denote it is in safety mode. Safety mode means that the scalpel is not fully extended. This feature avoids injuries when the blade is half way extended as it retracts back and does not pierce through. It has to be fully extended to use.
- **Safety Lock:** if you completely retract back it makes a loud click denoting it is in safety lock position. Safety lock position means that the blade cannot be reused once it is in lock position. This feature forces the doctor to dispose the product and minimizes the risk of transmissibility of blood-borne infectious agents.
- Compatible for a left handed as well as right handed surgeon.
- Safety feature in pre use and post use of product.
- One-handed activation allowing the clinician to focus on the surgical site.
- The retractable design protects the blade from exposure, minimizing the risk of accidental cuts.
- Universal size handle that fits all types of hand sizes, and ideal for kits, trays & packs.
- One handle size fits all blade sizes.

Details:

Product Name:	Universal Safety Scalpel
Size of Inner Box (cm):	17 x 8 x 3.5
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	42 x 35.5 x 19.5
Per 20 feet container:	280 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
9	PSF9	3	18	PSF18	4
10	PSF10	3	19	PSF19	4
10A	PSF10A	3	20	PSF20	4
11	PSF11	3	21	PSF21	4
11K	PSF11K	3	22	PSF22	4
13	PSF13	3	22A	PSF22A	4
14	PSF14	3	23	PSF 23	4
15	PSF15	3	24	PSF24	4
15B	PSF15B	3	25	PSF25	4
15C	PSF15C	3	MSC	PSFMSC	3
16	PSF16	3			
17	PSF17	3			



Biopsy Punch

Stainless Steel blade

Description

Biopsy punch is a sharp rounded blade tip used to remove tissue to be tested during procedures. The instrument is sterile with a round plastic handle. The punch is available in a variety of tip sizes and is sold in both sterile single pack as well as in bulk quantities.


Finely ground seamless, stainless steel blade ensures safe, clean and easy sample extraction. The sleek ribbed handle clear view of the biopsy site, comfort and control.

Virgin grade ABS plastic provides strength and sturdiness.

Details:

Product Name:	Biopsy Punch
Size of Inner Box (cm):	17 x 8 x 3.5
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	42 x 35.5 x 19.5
Per 20 feet container:	635 Carton

Size mm	Product Code	Product Code
1.0	BPS1.0	BPP1.0
1.5	BPS1.5	BPP1.5
2.0	BPS2.0	BPP2.0
2.5	BPS2.5	BPP2.5
3.0	BPS3.0	BPP3.0
3.5	BPS3.5	BPP3.5
4.0	BPS4.0	BPP4.0
5.0	BPS5.0	BPP5.0
6.0	BPS6.0	BPP6.0
7.0	BPS7.0	BPP7.0
8.0	BPS8.0	BPP8.0



BPS

BPS Features

- Ergonomically designed for precise skin excisions
- Unique, ribbed, easy-grip barrel and pressure assistance platform
- Sterile and ready to use.
- Wide range of sizes



BPP

BPP Features

- Premium quality punch
- Hard blister packaging
- Option to print custom brand on handle
- Pressure assistance platform and ribbed handle
- Micro grinding finish



Dermal Curette

Stainless Steel blade

Description

The curette has a firm place in the range of dermatological, diagnostic and therapeutic instruments. It is one of the most versatile instruments. For that reason, it is increasingly also used in other disciplines, such as plastic surgery.


The aim of a curette is the removal of a superficial skin changes with minimal damage to the surrounding tissue.

The curette is extremely sharp thus ensuring precise and gentle tissue removal. The outstanding cutting ability makes working with it much easier compared with conventional products.

Details:

Product Name:	Dermal Curette
Size of Inner Box (cm):	21 x 7 x 3
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	44 x 37 x 17
Per 20 feet container:	475 Carton

Size mm	Product Code	Product Code	Product Code	Product Code
2.0	DCS2.0	DCP2.0	DCK2.0	DDC2.0
3.0	DCS3.0	DCP3.0	DCK3.0	DDC3.0
4.0	DCS4.0	DCP4.0	DCK4.0	DDC4.0
5.0	DCS5.0	DCP5.0	DCK5.0	DDC5.0
7.0	DCS7.0	DCP7.0	DCK7.0	DDC7.0



DCS Features

- Seamless stainless steel sharp cutting edge
- Easy dynamic gripping handle
- Sterile and ready to use
- Sharp side indicator
- Wide range of sizes



DCK Features

- Sturdy & Durable Stainless Steel
- Symmetrically Ground Beveled Edge
- Ultra Sharp Inside Grinding
- Sharp side indicator
- Ribbed on all sides
- Triangular Handle
- Sterile and ready to use
- Wide range of sizes




DCP Features

- Long 'Pen Type' Grip
- Branding Labeling area
- Ribbed handle for easy control
- Minimum damage to surrounding tissues
- Sharp side indicator
- Sterile and ready to use
- Wide range of sizes



DDC Features

- Seamless stainless steel sharp cutting edge
- One in metal with a sharp edge
- Sharp side indicator
- Covered by a protective cap
- Ultra Sharp Inside Grinding
- Sterile and ready to use
- Wide range of sizes



Ophthalmic Knives

Round-Wire Blades



Slit-Blades (Keratome) For Phaco Stab Incision

Product No. Description

P-912301 2.3 MM Sharp Tip Straight
P-912501 2.5 MM Sharp Tip Straiaht
P-912601 2.65 MM Sharp Tip Straight
P-912801 2.8 MM Sharp Tip Straiaht
P-912901 3.0 MM Sharp Tip Straight
P-913201 3.2 MM Sharp Tip Straiaht
P-913501 3.5 MM Sharp Tip Straight
P-915101 5.1 MM Sharp Tip Straiaht

Product No. Description

P-912361 2.3 MM Sharp Tip Angled 45° 19 Gauge
P-912561 2.5 MM Sharp Tip Angled 45° 19 Gauge
P-912661 2.65 MM Sharp Tip Angled 45° 19 Gauge
P-912861 2.8 MM Sharp Tip Angled 45° 19 Gauge
P-912961 3.0 MM Sharp Tip Angled 45° 19 Gauge
P-913261 3.2 MM Sharp Tip Angled 45° 19 Gauge
P-913561 3.5 MM Sharp Tip Angled 45° 19 Gauge
P-915061 5.1 MM Sharp Tip Angled 45° 19 Gauge

Double Bevel

P-912808 2.8 MM Sharp Tip Straight
P-912908 3.0 MM Sharp Tip Straight
P-913208 3.2 MM Sharp Tip Straight

Double Bevel

P-912868 2.8 MM Sharp Tip Angled 45° 19 Gauge
P-912968 3.0 MM Sharp Tip Angled 45° 19 Gauge
P-913268 3.2 MM Sharp Tip Angled 45° 19 Gauge



Slit-Blades (N-Larger) For I.O.L. Enlarsins

Product No. Description

P-914001 4.0 MM Blunt Tip Straight
P-915201 5.2 MM Blunt Tip Straight
P-915501 5.5 MM Blunt Tip Straight
P-916001 6.0 MM Blunt Tip Straight
P-916201 6.2 MM Blunt Tip Straight

Product No. Description

P-914061 4.0 MM Blunt Tip Angled 45° 16 Gauge
P-915561 5.5 MM Blunt Tip Angled 45° 15 Gauge
P-916061 6.0 MM Blunt Tip Angled 45° 15 Gauge
P-916261 6.2 MM Blunt Tip Angled 45° 15 Gauge



Crescent Blades For Groove Or Pocket Or Tunnel Incision

Product No. Description

P-950001 Crescent Sharp Tip Straight

Product No. Description

P-950002 Crescent Bevel Up Angled 45° 20 Gauge
P-950003 Crescent Bevel Down Angled 45° 20 Gauge

Double Bevel

P-950004 Crescent Angled 45° 20 Gauge
P-950005 Sharp Tip Straight



Lance Tip Blades For Initial Incision

Product No. Description

P-931501 15° Degree 24 Gauge
P-933001 30° Degree 24 Gauge
P-934501 45° Degree 24 Gauge

Round-Wire Blades



MVR Blades For Water Tight Self Sealing Incision

Product No. Description

P-975559 MVR 24 Gauge
P-975560 MVR 20 Gauge
P-975561 MVR 19 Gauge

Product No. Description

AC STILETIO BLADES
For Blumenthal or Side Port Incision
P-985560 MVR 20 Gauge Angled 45°
P-985561 MVR 19 Gauge Angled 45°



Miniature Blades Alternative To General Surgery For Precise Incision

Product No. Description

SPOON BLADES
For Groove or Pocket or Tunnel Incision

Product No. Description

P-6821 3.0 MM Angled 60°
P-6821E 3.0 MM Extended Neck Angled 60°



Miniature Blades Alternative To General Surgery For Precise Incision

Product No. Description

SCLERAL BLADES
For Limbal and Scleral Incision
P-5700 Wide Tip Blade
P-5710 Narrow Tip Blade



Details:

Product Name:	Ophthalmic Knives
Size of Inner Box (cm):	20 x 5 x 4
Unit/Box:	6 Pcs
Export carton:	25 Boxes x 2 Outer
Export Carton size (cm):	42.5 x 27.5 x 22.5
Per 20 feet container:	1000 Carton

Chisel Micro Surgical Blades

Fine Stainless Steel

Description

FINE (CHISEL) Surgical Blades Stainless Steel are available in the following sizes: Nos. Each blade is protected by VCI paper in a peel-open foil sachet. Sterilized by Gamma Radiation of 2.5 M. Rads.

Details:

Product Name:	Chisel Blades (Fine)
Size of Inner Box (cm):	8.2 x 3.2 x 2.2
Unit/Box:	25 Pcs
Outer Boxes:	25 Boxes x 25 Units
Export carton:	25 Boxes x 10 Outer
Export Carton size (cm):	47.5 x 39 x 17
Per 20 feet container:	825 Carton

Size	Product Code	Handle
61	PSM61	PSMH
61V	PSM61V	PSMH
62	PSM62	PSMH
62V	PSM62V	PSMH
63	PSM63	PSMH
64	PSM64	PSMH
65	PSM65	PSMH
66	PSM66	PSMH
67	PSM67	PSMH
68	PSM68	PSMH
69	PSM69	PSMH
69B	PSM69B	PSMH
90	PSM90	PSMH
91	PSM91	PSMH



Skin Graft Blade

Carbon Steel

Description

Skin Graft Blade, Simplex, Duplex


Surgical Skin Graft Blades are 158 mm long, individually wrapped in VCI paper in single peel packets in packs of 10.

These blades can be used with Watson, Braithwaite or Cobbett Knives.

Details:

Product Name:	Skin Graft Blades
Size of Inner Box (cm):	20.7 x 3.6 x 2
Unit/Box:	10 Pcs
Outer Boxes:	10 Boxes x 10 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	44 x 24.5 x 22.5
Per 20 feet container:	1100 Carton

Size	Carbon Steel	Stainless Steel
SG	PCSG	PSSG
SGS	PCSGS	PSSGD
SGD	PCSGS	PSSGD




Blood Lancet

Description

Supplied sterile and individually packed. Lancet comes in a pack of 200 pieces per box. the tip is specially designed for easy pricking, which minimizes patient trauma.

Details:

Product Name:	Blood Lancet
Product Code:	PSSBL
Size of Inner Box (cm.):	12.5 x 7 x 20
Unit/Box:	200
Export carton:	100 Boxes x 200 Units
Export Carton size (cm.):	38.5 x 26.5 x 15.5
Per 20 feet container:	1800 Cartons



Myringotomy Knife Blades

Description

Myringotomy is a surgical procedure of the eardrum or tympanic membrane. The procedure is performed by making a small incision with a myringotomy knife through the layers of tympanic membrane.

Features:

- Optimised micrometre-accurate precision grinding.
- Double edged allowing for cutting movements in both directions.
- Ultra-sharp cutting edges of superb durability.
- Consistently sharp, angled blade facilitates a clean cut.
- Spear tip, Narrow shaft.
- Lance tip, Downward cut, Narrow shaft .
- With and Without Handle



Details:

Product Name:	Myringotomy Knife Blade
Size of Inner Box (cm):	24.2 x 8 x 2.5
Unit/Box:	5 Pcs
Outer Boxes:	--
Export carton:	50 Boxes
Export Carton size (cm):	50.5 x 42 x 14
Per 20 feet container:	1000 Carton

Size	Product Code
Lance	PSSML
Spear	PSSMS



Face Mask - 3 Ply/ Type IIR



3 PLY Mask - Non Surgical



Smooth Fiber
Material



SMS Fabric

3 PLY Mask - Surgical



Meltblown
Fabric



UV Cleaned

FEATURES




Adjustable Nose
Clip



Ultrasonic
Sealing



One Size
Fits Most



Seamless
Crimping



Comfortable
& Breathable



Quality
Material



3 Layer
Protection



Anti Droplets

Details:

Product Name:	3 PLY/ Type IIR Mask - Non Surgical
Product Code	P03N
Size of Inner Box (cm.):	23 x 11 x 19.5
Unit/Box:	100
Export Carton size (cm.):	57 x 48 x 41
Per 20 feet container:	230 Carton

Details:

Product Name:	3 PLY/ Type IIR Mask - Surgical
Product Code	P03M
Size of Inner Box (cm.):	23 x 10.5 x 10
Unit/Box:	50
Export Carton size (cm.):	55 x 48 x 42
Per 20 feet container:	230 Carton

Face Mask - N95/FFP2



without Respirator

with earloop



UV Cleaned



without Respirator

with headloop



Filters Particulates



BFE≥95%



with Respirator

with earloop



PFE≥95%

FEATURES



Adjustable Noise Clip



Soft Elastic Ear Loops



High Quality Breathable
Material and Valve



Seamless Flow of
Exhaled Air

1st Layer

Fine mesh.
Filters unwanted particulates

3rd Layer

Hot air mesh.
Blocks smaller particulates

2nd Layer

Heavy meltblown filter mesh.
Blocks odour, virus, pollution

4th Layer

Extra meltblown.
Filters remaining particulates

5th Layer

Extra soft non woven
layer to fit the face

Details:

Product Name:	N95/FFP2 without Respirator
Product Code:	P95W
Size of Inner Box (cm.):	30.5 x 13.5 x 17.5
Unit/Box:	50
Export Carton size (cm.):	63.5 x 56 x 55
Per 20 feet container:	140 Carton

Details:

Product Name:	N95/FFP2 with Respirator
Product Code:	P95R
Size of Inner Box (cm.):	24.5 x 13.5 x 17.5
Unit/Box:	25
Export Carton size (cm.):	51 x 42 x 54.5
Per 20 feet container:	220 Carton

Disposable Adult Diapers


For Urinary Incontinence Management

Features:

- Dry-Fast Core, Super absorbent.
- Moisture repellent.
- Security barriers.
- Waist Elastication and curved leg elastic.
- Superfit Band Wetness Indicator.
- Three absorbency levels and three sizes.
- Compact and convenient packaging.
- Also available in options of double cotton core & ADL.

Benefits:

- Faster absorption, drier surface.
- Better containment of faces.
- Reduces urine leakage.
- Improves comfort.
- Comfortable fit.
- Easier to fix and adjust tapes.
- A Guide for carers.
- Meets individual requirements with cost effective care.
- Easy to handle less storage space required.



Details:

Product Name:	Adult Diapers (Medium)
Product Code:	AKMP 04
Size of Poly Bag (cm.):	20 x 23.0 x 17.5
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size (cm.):	46 x 35.5 x 60
Per 20 feet container:	260 Cartons

Details:

Product Name:	Adult Diapers (Large)
Product Code:	AKLP 06L
Size of Poly Bag (cm.):	22 x 24 x 17.5
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size (cm.):	48.5 x 35.5 x 66
Per 20 feet container:	225 Cartons

Details:

Product Name:	Adult Diapers (Extra Large)
Product Code:	AKXLP 07
Size of Poly Bag (cm.):	23.5 x 24 x 17.5
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size (cm.):	49 x 35.5 x 70
Per 20 feet container:	205 Cartons



Under Pads

AroKleen

Description

- **Anti-Leakage Design:**
The product has extra width and length designed for better absorption.
- **Dry And Comfortable:**
It has large absorbing padding using the reticular embossed design to introduce the water diversion towards every side keep the surface comfortable.
- **Skin Friendly:**
Using the super-nonwoven as the top layer, it feels comfortable and prevents skin rashes.
- **International Quality:**
Strict controls of product material and production assuring safety and health.

Details:

Product Code:	ABM60x60
Size of Poly Bag (cm.):	29 x 20 x 10
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size (cm.):	41 x 31 x 59
Per 20 feet container:	350 Cartons


Details:

Product Code:	US90x60
Size of Poly Bag (cm.):	29 x 20 x 10
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size:	41 x 31 x 59
Per 20 feet container:	350 Cartons

Details:

Product Code:	PP30P
Size of Poly Bag (cm.):	26.5 x 18 x 27
Unit/Box:	30
Export carton:	12 Bags x 30 Units
Export Carton size (cm.):	60 x 36 x 29
Per 20 feet container:	400 Cartons

5 LAYERED LEAK PROOF TECHNOLOGY!



The Advantages of Using PupPee Pads

Our pads provide worry-free training for your puppy & peace of mind for you. They also can be used as an alternative to outdoor "relief" for your adult dog.

- Locks in wetness with quilted top layer
- Ensures maximum absorption power with 2 moisture-locking layers
- Protects floor & carpets with leak-proof backing
- Provides easy training with built-in attractant
- Protect your floors
- Easy cleanup & disposal
- Quilted for extra absorbency
- Less accidents means less disciplining
- Dog attractant reinforces training
- Durable, tear resistant



Disposable Adult Diapers

Premium



ADL (Acquisition Distribution Layer)



Details:

Product Name:	APM10
Size of Poly Bag (cm.):	15X18.5X18.5
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size (cm.):	47.5X42X38.5
Per 20 feet container:	330 Cartons

Economy



Details:

Product Code:	APL11
Size of Poly Bag (cm.):	16X19X21
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size:	47.5X42X38.5
Per 20 feet container:	330 Cartons

Details:

Product Name:	APX12
Size of Poly Bag (cm.):	20X16X21.5
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size (cm.):	47.5X42X38.5
Per 20 feet container:	330 Cartons



Disposable Prep Razor

Paraprep

Features

- Safe contoured handle with sure grip ridges provide precise fingertips control.
- Easy to handle lies in front of blades thus requiring only light gentle pull instead of pressure.
- Coated stainless steel blades shaves close without nick or irritation.
- Open design provides optimal view of shaving area.
- Compact and convenient packaging.
- Cost effective, no need for more than one razor per procedure.



Details:

Product Code:	P/DR20210320
Size of Poly Bag (Cm.):	9.5x5.7
Inner Box:	20x12x6.5
Unit/Box:	50
Export carton:	12 Bags x 50 Units
Export Carton size:	41.5X27.5X35
Per 20 feet container:	650 Cartons



with Two Side Blades



For All Body Areas



Smooth Hair Removal



100% Shower Proof



Global Standard Manufacturing



Unique Skin Protection



Stainless Steel Blade

Details:

Product Code:	P/DR20211020
Size of Poly Bag (Cm.):	8.5x15.5
Inner Box:	28.5x19X26
Unit/Box:	100
Export carton:	12 Bags x 100 Units
Export Carton size:	41.5X27.5X35
Per 20 feet container:	650 Cartons



Lung Exerciser

Spirometer

Features:

- Three color balls for easier identification of inhale capacity
- Improve maximal exercise capacity
- Made from transparent material
- Compact, convenient and safe



Details:

Product Code:	PEL 01
Size of Box (Cm.):	14.2x7x15.7
Unit/Box:	01
Export carton:	12 Bags x 01 Units

Pulse Oximeter

Features

- 2pcs standard AAA1.5V alkaline battery supporting 20 hours continuous operation.
- Applicable to Hospital, Clinic, Family, Oxygen Bar, Sporting center etc.
- True color high resolution OLED/TFT LCD/LED.
- Light in weight, simple in operation, easy to use & carry.
- Adjust the parameters in friendly menu.
- Soft anti-ambient light interference.
- Low power consumption.
- Optional Pouch.



Details:

Product Code:	PPOM-09
Size of Box (Cm.):	8.5x5x6
Unit/Box:	01
Export carton:	12 Bags x 01 Units

Infrared Forehead Thermometer

Features

- Non-contact Forehead
- LCD Backlight
- C/F Unit Switchable



16 Memories



High Temperature
Alarm



1s Fast
Reading



Details:

Product Code:	PPOM-09
Size of Box (Cm.):	10.5X4.65X16.5
Unit/Box:	01
Export carton:	12 Bags x 01 Units

Bed Bath Wipes

AroKleen

Description

- AroKleen Bed Bath Patient Wipes are Paraben Free, Alcohol Free, SLS Free, Mircoavable.
- The fast acting formula, with a two minute kill time loosens and suspends blood and bodily fluids, non-organic soil and other cellular debris for easy cleaning.

CLEANS, DISINFECTS & DEODORIZES

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



Skin Friendly
pH



Aloe Vera



Enriched
With Glycerin

USAGES



Cleansing



Refreshing



Moisturizing



Hospital Care

Details:

Product Code:	BP2202
Size of Poly Bag (Cm.):	10X8.5X24
Unit/Box:	10
Export carton:	12 Bags x 10 Units
Export Carton size:	43 x 24 x 23
Per 20 feet container:	1100 Cartons



Baby Wipes

AroKleen

Description

- AroKleen Baby Wipes are large, durable, disposable clothes that are alcohol and latex free. Our Aloe Vera Extract, Vitamin E lotion-based formula that gently cleans and keeps babies' skin smooth and moisturized.

SPECIALLY FOR BABIES

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



No Added
Parabens



No Added
Alcohol



Extra Thick &
Extra Soft



Silicon Free



Enriched With
Vitamin E



pH
Balanced



Aloe Vera
Extract



Details:

Product Code:	BB23009
Size of Poly Bag (Cm.):	10X10.3X24.5
Unit/Box:	72
Export carton:	12 Bags x 72Units
Export Carton size:	36 x 32 x 23.5
Per 20 feet container:	950 Cartons



Cleansing & Refreshing Wipes

SPECIALLY FOR CLEANSING

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



No Added
Alcohol



Skin Friendly
pH



Enriched With
Vitamin C

Details:

Product Code:	BN2202
Size of Poly Bag (Cm.):	10X9.2X23.5
Unit/Box:	30
Export carton:	12 Bags x 30 Units
Export Carton size:	43 x 24 x 23
Per 20 feet container:	1100 Cartons



Germ Protection Wipes

AroKleen

Description

- AroKleen Germ Protection Wipes are more effective and cost less than sprays and towels, making them the right wipe to use in any location. Every area where staff, visitors, guests, customers and the general public convene, there are surfaces that need to be regularly cleaned and freed of dirt, grime, harmful germs and bacteria.

PROTECTION FROM GERMS

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



Skin Friendly
pH



Aloe Vera

USAGES



Hands



Refreshing



Restaurant



Car Interior



Device



Toys

Details:

Product Code:	PGP01
Size of Poly Bag (Cm.):	10X8.5X24
Unit/Box:	48
Export carton:	12 Bags x 48 Units
Export Carton size:	37 x 32 x 22.5
Per 20 feet container:	950 Cartons



Pet Wipes

SPECIALLY FOR CATS & DOGS

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



Powerful Odour



pH Balanced



Olive Oil



Lavender
Essential Oil



Enriched With
Vitamin E



Cruelty Free



Vet & Groomer
Recommended

Details:

Product Code:	PGP01
Size of Poly Bag (Cm.):	10X10X24
Unit/Box:	48
Export carton:	12 Bags x 48 Units
Export Carton size:	37 x 32 x 22.5
Per 20 feet container:	950 Cartons



Hot Water Bottle



Joint Pain



Sports Pain



Arthritis Pain



Rheumatic Pain



Back Aches



Stomach Aches



Muscular Cramps



Muscle Pull



Regular Hot Water
Bottle



BS Standard Hot
Water Bottle

DNV-GL

EC CERTIFICATE

Full Quality Assurance System

Certificate No.: 10000400574 PA-NA-IND Project No.: PRC-213992-2010-PRC-IND Valid Until: 27-05-2024

This is to certify that the quality system of:

Paramount Surgimed Ltd.

Works: **A-106, RIICO Industrial Area, Bhiwadi – 301 019, District Alwar, Rajasthan, India**

For design, production and final product inspection/testing of:

STERILE SURGICAL DISPOSABLE MEDICAL DEVICES

Has been assessed with respect to:

THE CONFORMITY ASSESSMENT PROCEDURE DESCRIBED IN ANNEX II EXCLUDING SECTION 4 OF COUNCIL DIRECTIVE 93/42/EEC ON MEDICAL DEVICES, AS AMENDED

and found to comply.

Further details of the product(s) and conditions for certification are given overleaf.

Place and date:
Havli, 15 September 2020

For:
DNV GL PRESAFE AS
Notified Body No.: 2460

Mariann Jeremiasen

The certificate is digitally verified by blockchain technology. For more info, see www.dnvgl.com/assessing/certificates-on-the-blockchain.html

Notice: The Certificate is subject to terms and conditions as set out in the Certification Agreement. Failure to comply may render this Certificate invalid.
NOTIFIED BODY: DNV GL PRESAFE AS, Verkearsveien 5, N-1363 Havli, Norway - Registered Enterprise No: NO 987 067 401 MVA.

DNV-GL-078-A Rev. 0.0 Page 1 of 4

EC Certificate

DNV-GL

Certificate No.: 10000400574 PA-NA-IND Project No.: PRC-213992-2010-PRC-IND Valid Until: 27-05-2024

Jurisdiction

Application of Council Directive 93/42/EEC of 14 June 1993, adopted as "Forskrift om Medisinsk Utstyr" by the Norwegian Ministry of Health and Care Services.

Certificate history:

Revision	Description	Issue Date
0.0	Original Certificate	2020-09-15

Products covered by this Certificate:

Product Description	Product Name	Class
Sterile Curette Dermal	Curette Dermal Size: 2, 3, 4, 5 & 7	Ila
Sterile Stitch cutters in Carbon Steel and Stainless Steel	Stitch Cutters Long, Short, Mini	Ila
Sterile Surgical Blades in Carbon Steel and Stainless Steel	Surgical blades 1, 2, 3, 4, 5, 6, 8, 9, 10, 10A, 11, 11K, 12, 13, 14, 15, 15B, 15C, 15D, 15S, 16, 17, 18, 19, 20, 21, 22, 22A, 23, 24, 25, 36, 40, 40B, 60, 60B, 1R, 2R, 3R, 1V, 2V, 3V, 11P, 12D, 24D, 34, 36D	Ila
Sterile Disposable Scalpels in Carbon Steel and Stainless Steel	Disposable Scalpels 1, 2, 3, 4, 5, 6, 8, 9, 10, 10A, 11, 11K, 12, 13, 14, 15, 15B, 15C, 15D, 15S, 16, 17, 18, 19, 20, 21, 22, 22A, 23, 24, 25, 36, 1R, 2R, 3R, 1V, 2V, 3V, 11P, 12D, 24D, 34, 36D	Ila
Sterile Safety Scalpels in Carbon Steel and Stainless Steel	Safety Scalpels 6, 9, 10, 10A, 11, 11K, 12, 13, 14, 15, 15B, 15C, 15D, 15S, 16, 17, 18, 19, 20, 21, 22, 22A, 23, 24, 25, 36, 1R, 2R, 3R, 1V, 2V, 3V, 11P, 12D, 24D, 34, 36D	Ila
Sterile Fine Blades / Chisel Blade / Microsurgery blades	61, 62, 63, 64, 65, 66, 67, 68, 69, 69B, 61V, 62V, 60, 91	Ila
Sterile Ophthalmic Blades	Keratome: P-912301, P-912501, P-912601, P-912801, P-912901, P-913201, P-913501, P-915001, P-912361, P-912561, P-912661, P-912861, P-912961, P-913261, P-913561, P-915061, P-912808, P-912908, P-913208, P-912868, P-912968, P-913268, P-914001, P-915201.	Ila

NOTIFIED BODY: DNV GL PRESAFE AS, Verkearsveien 5, N-1363 Havli, Norway - Registered Enterprise No: NO 987 067 401 MVA.

DNV-GL-078-A Rev. 0.0 Page 2 of 4

EC Certificate

DNV-GL

Certificate No.: 10000400574 PA-NA-IND Project No.: PRC-213992-2010-PRC-IND Valid Until: 27-05-2024

	P-915501, P-916001, P-916201, P-914061, P-915161, P-915561, P-916061, P-916261 Crescent: P-950001, P-950002, P-950003, P-950004, P-950005 Lance Tip: P-931501, P-933001, P-934501, P-913501, P-915101, P-915161, P-914101, P-914161, P-915261 MVR: P-975559, P-975560, P-975561, P-985560, P-985561 Spoon: P-6820, P-6821, P-6821E Scleral: P-5700, P-5710	
Sterile Biopsy Punches	Biopsy Punches: 1mm, 1.5mm, 2mm, 2.5mm, 3mm, 3.5mm, 4mm, 5mm, 6mm, 7mm, 8mm, 10mm, 12mm, 15mm	Ila
Sterile Skin Graft Blades	Simplex, Duplex	Ila
Sterile Blood Lancet	Standard	Ila
Sterile Myringotomy	Lance and Spear	Ila

The complete list of devices is filed with the Notified Body

Sites covered by this certificate

Site Name	Address
Paramount Surgimed Ltd.	A-106, RIICO Industrial Area, Bhiwadi – 301 019, District Alwar, Rajasthan, India

EU Representative

Medical Device Safety Service GmbH,
Schiffgraben 41, D-30175, Hannover, Germany.

NOTIFIED BODY: DNV GL PRESAFE AS, Verkearsveien 5, N-1363 Havli, Norway - Registered Enterprise No: NO 987 067 401 MVA.

DNV-GL-078-A Rev. 0.0 Page 3 of 4

EC Certificate

DNV-GL

Certificate No.: 10000400574 PA-NA-IND Project No.: PRC-213992-2010-PRC-IND Valid Until: 27-05-2024

Terms and conditions

The certificate is subject to the following terms and conditions:

- Any producer (see 2001/95/EC for a precise definition) is liable for damage caused by a defect in his product(s), in accordance with directive 85/374/EEC, as amended, concerning liability of defective products.
- The certificate is only valid for the products and/or manufacturing premises listed above.
- The Manufacturer shall fulfil the obligations arising out of the quality system as approved and uphold it so that it remains adequate and efficient.
- The Manufacturer shall inform Presafe of any intended updating of the quality system and Presafe will assess the changes and decide if the certificate remains valid.
- Periodical audits will be held, in order to verify that the Manufacturer maintains and applies the quality system. Presafe reserves the right, on a spot basis or based on suspicion, to pay unannounced visits.

The following may render this Certificate invalid:

- Changes in the quality system affecting production.
- Periodical audits not held within the allowed time window.

Conformity declaration and marking of product

When meeting with the terms and conditions above, the producer may draw up an EC declaration of conformity and legally affix the CE mark followed by the Notified Body identification number of Presafe.

End of Certificate

NOTIFIED BODY: DNV GL PRESAFE AS, Verkearsveien 5, N-1363 Havli, Norway - Registered Enterprise No: NO 987 067 401 MVA.

DNV-GL-078-A Rev. 0.0 Page 4 of 4

EC Certificate



Certificate no.:	Initial certification date:	Valid:
347843.2017.AQ.IND.NA.PS	03 March 2018	04 March 2022 – 03 March 2025

has been found to conform to the Quality Management System standard:
ISO 13485:2016 / EN ISO 13485:2016

Place and date:
Month: 07 March 2020



For the issuing office:
DNV Product Assurance AS
 Veritetsveien 2, 1163 Hovik, Norway

Cecilie Gudesan Torp
 Associate Professor, Department of Health, Behavior, and Society, Johns Hopkins University

Lack of fulfillment of conditions as set out in the Certification Agreement may render this Certificate invalid.

ISO : 13485



Certificate no.: 247943-2017-AQ-IND-NA-PS
Place and date: Hanoi, 07 March 2022

Paramount Surgimed Ltd.
Locations included in the certification are as follows:

Site Name	Site Address	Site Scope
Paramount Surgimed Ltd.	Plot - 1, L.S.C., Okhla Industrial Area, Phase - II, New Delhi - 110 025, India	Design & Development, Sales & marketing of Sterile / Unsterile Surgical Blades, Disposable Scalpels, Disposable Retractable Safety Scalpels, Stitch Cutters, Skin Graft Blades in Carbon Steel and Stainless Steel and Blood Lancets, Biopsy Punches, Ophthalmic Blades, Fire Blades/Chisel Blade/Microsurgery Blades & Curette Dermal in Stainless Steel
Paramount Surgimed Ltd.	Site 2 - A-106, RICO Industrial Area, Bhivadi, District Awar - 301 019, Rajasthan, India	Manufacture and Supply of Sterile / Unsterile Surgical Blades, Disposable Scalpels, Disposable Retractable Safety Scalpels, Stitch Cutters, Skin Graft Blades in Carbon Steel and Stainless Steel and Blood Lancets, Biopsy Punches, Ophthalmic Blades, Fire Blades/Chisel Blade/Microsurgery Blades & Curette Dermal in Stainless Steel
Paramount Surgimed Ltd.	Site 3 - Plot No. 33, Sector-14, Phase-2, MT, Bawal - 123501, (Haryana), India	Manufacture and Supply of Sterile / Unsterile Surgical Blades, Disposable Scalpels, Disposable Retractable Safety Scalpels, Stitch Cutters, Skin Graft Blades in Carbon Steel and Stainless Steel and Blood Lancets, Biopsy Punches, Ophthalmic Blades, Fire Blades/Chisel Blade/Microsurgery Blades & Curette Dermal in Stainless Steel

Lack of fulfillment of conditions as set out in the Certification Agreement may render this Certificate invalid.

Page 2 of 2

ISO : 13485



This is to Certify that

H.Q.: 1, L.S.C., Okhla Industrial Area, Phase II, New Delhi - 110 020, India.
Works: A - 106, RIICO Industrial Area, Bhiwadi, District Alwar - 301019 (Rajasthan), India.
Works: Plot No. 33, Sector - 14, Phase - 2, IMT Bawal - 123501 (Haryana), India.

has been assessed by RSRV and found to comply with the requirements of

For the following Scope:

Design & Development, Manufacturing, Marketing and Supply of Medical Products like Adult Diapers, Disposable Undershorts, Baby Diapers, Sanitary Napkins, Under Pads, 2 Ply Face Masks, 3 Ply Face Masks, KN95 Face Masks, N95 Face Masks/ FFP2/FF3, PPE Kit, Prep Razors, Baby Wipes, Bed Bath Wipes, Cleansing & Refreshing Wipes, Germ Protection Wipes, Disinfectant Wipes, Makeup Removal Wipes, Facial Cleansing Wipes, Dermal Blades and Suction Needle.

Certificate No. : RSBV-2386-6326

Date of This Certificate Issue: 03/06/2023
1st Annual surveillance audit due on: 02/06/2024
2nd Annual surveillance audit due on: 02/06/2025
Validatum / Date of Expiry: 02/06/2026

Royal Stancert B.V.

Rept. ENR06 : Wageningen UR, 1014 AH Alameda, The Netherlands.

This certificate remains the property of **Braskem S.A.** and must be returned whenever demanded the validity of this certificate can be verified at <http://www.insitesteel.org>. **Braskem S.A.** is an independent system, product and personnel assessment body accredited by Global Trust Accreditation (Atlanta, Georgia, U.S.A. - 014). Email: info@insitesteel.org



mno
LABORATUVARLAR

Belge No / Certificate No	19206266
Belge/Onaylama Tarihi - Bir Sanyaklı Belge Tarihi / Certification Date / Certificate Validity Date	10.07.2023/07.07.2028
Belge Geçerlilik Tarihi / Document Validity Period	15 Yıl / 5 Years
Firma Ünvanı ve Adresi / Company Name and Address	Paramount Servisler Ltd

Marka / Model / Brand / Model
 Deskripsi / Description
 Modifikasi/Warranty / Modification / Warranty

Tehnik Egeyendilme Hapur No 1

Form Type / Product Type	Form Type / Product Type
--------------------------	--------------------------

• EN 149:2001 - A1 2009 Schmanke diaphragm filter - Pasopokura karp korima anagh filtri yama sandakali/ Respiratory protective device - Filtering half mask to protect against particles

Çizimdeki Malzeme Bilgisi / Product Material Information: Aşılıcam model tiruloru kamay, elasto kayış, boros lifleri ve diğer kısımları kollektörler tarafından üretilmiştir / Aşılıcam model products are manufactured using fabric, elastic straps, boros fibers and other parts by collectors

Yoklav, ARJN
01/07/2022
Karat Field / Agronomer

Class A&Z,
8/10/2023
Sivak Muthu / General manager



Adres: Alameda da Universidade, Caixa Postal 12461, Curitiba, PR, Brasil
Tel: 011 3361-0100 Fax: 011 3361-0101 www.ufrpr.br



GMP

CE Mask



1, L.S.C., Okhla Industrial Area, Phase-II, New Delhi-110 020, India

Tel +91-11-46436600, 46436601, 46436666, 26389812, 13

Fax +91-11-26389815, 41616555

sales@paramountblades.com

export2@paramountblades.com


www.paramountsurgimed.com

(Disclaimer- All measurements in weight & size are indicative for reference only.)

PSL/CTLG/01, Issue No. 01, Issue Date : 26/11/2020, Rev. No. 03, Rev. Date : 29/04/2023

topedge

PRODUCT CATALOGUE





topedge

Surgical Blades	02
-----------------	----

Stitch Cutters	04
----------------	----

Chisel Blades	06
---------------	----

Gouge Blades & Handles	07
------------------------	----

Disposable Scalpels	08
---------------------	----

Universal Safety Scalpel	09
--------------------------	----

Biopsy Punch	10
--------------	----

Dermal Curette	11
----------------	----



Ophthalmic Knives	12
-------------------	----

Surgical Blades

Carbon Steel & Stainless Steel

Features

- Stainless Steel Surgical Blades and Carbon Steel Surgical Blades
- Packed in peel open pouches with a VCI Protective inner line
- Stainless Steel (Grade F as Per ISO 7153-1:2016)
- Carbon Steel (Grade - As Per BS 2982:1992)
- Complies the Hardness requirement of BS 2982:1992 standard
- Complies the fitment requirement of ISO 7740:1985 standard requirement.
- Shelf life of 5 years



Details :

Product Name :	Surgical Blades
Size of Inner Box (cm) :	8.5 x 8.1 x 2.7
Unit/Box :	100 pcs
Outer Boxes :	10 Boxes x 100 Units
Export carton :	10 Boxes x 10 Outer
Export Carton size (cm) :	47 x 39 x 17
Per 20 feet container :	860 Carton

Post Mortem Blades

Size 60

Features

- Used in the fields of Pathology and Histopathology for autopsy
- It is used with handle No.8
- Available in both carbon and stainless steel
- Box of 50
- Individual peel pouch



Size 40

Features

- The PM40 range of blades are meant to be used with screw type heavy duty handles
- Heavy duty blade
- Box of 10



Size	Carbon Steel	Stainless Steel	Fitment No
9	TCS9	TSS9	3
10	TCS10	TSS10	3
10A	TCS10A	TSS10A	3
11	TCS11	TSS11	3
11K	TCS11K	TSS11K	3
11P	TCS11P	TSS11P	3
12	TCS12	TSS12	3
12D	TCS12D	TSS12D	3
13	TCS13	TSS13	3
14	TCS 14	TSS14	3
15	TCS15	TSS15	3
15B	TCS15B	TSS15B	3
15C	TCS15C	TSS15C	3
15D	TCS15D	TSS15D	3
15S	TCS15S	TSS15S	3
16	TCS16	TSS16	3
17	TCS17	TSS17	3
18	TCS18	TSS18	4
19	TCS19	TSS19	4
20	TCS20	TSS20	4
21	TCS21	TSS21	4
22	TCS22	TSS22	4
22A	TCS22A	TSS22A	4
23	TCS23	TSS23	4
24	TCS24	TSS24	4
24D	TCS24D	TSS24D	4
25	TCS25	TSS25	4
34	TCS34	TSS34	4
36	TCS36	TSS36	4
36D	TCS36D	TSS36D	4
40	TCS40	TSS40	-
40B	TCS40B	TSS40B	-
60	TCS60	TSS60	8
60B	TCS60B	TSS60B	8

Stitch Cutters

Carbon Steel & Stainless Steel

Description

Stitch Cutter is used for general suture removal. Appearing similar to a scalpel, this instrument requires no handle and is a simple and convenient, way to remove basic sutures. The slightly curved, flat, bull nose tip slides under the suture and when rotated, gently removes the suture from the site. Long stitch cutter is made of stainless steel and short stitch cutter is made of carbon steel. The mini stainless steel stitch cutter can be used with handle No.3L, 7 & 9. Alternatively, it is also available on pre-mounted ABS disposable plastic handle.

Features

- Five Year Shelf Life – Sterile
- Each blade is protected by VCI paper in a peel-open foil sachet.
- Convenience of single-use disposables.
- Sizes clearly marked on each blade.
- Each box contains 100 pcs

Size	Carbon Steel	Stainless Steel	Handle No.	Handle No.
LSC	--	TLSCSS	TSS3L	TSS9
SSC	TSSCCS	--	TSS3L	TSS9
MSC	--	TMSCSS	TSS3L	TSS7,9



Details :

Product Name :	Mini Stitch Cutter
Size of Inner Box (cm) :	8.5 x 8.1 x 2.7
Unit/Box :	100 Pcs
Outer Boxes :	10 Boxes x 100 Units
Export carton :	10 Boxes x 10 Outer
Export Carton size (cm) :	47 x 39 x 17
Per 20 feet container :	860 Carton

Details :

Product Name :	Short Stitch Cutter
Size of Inner Box (cm) :	11.1 x 8 x 3.5
Unit/Box :	100 Pcs
Outer Boxes :	10 Boxes x 100 Units
Export carton :	10 Boxes x 10 Outer
Export Carton size (cm) :	46.5 x 39.8 x 25
Per 20 feet container :	570 Carton



Details :

Product Name :	Long Stitch Cutter
Size of Inner Box (cm) :	14.8 x 7 x 3.1
Unit/Box :	100 Pcs
Outer Boxes :	10 Boxes x 100 Units
Export carton :	10 Boxes x 10 Outer
Export Carton size (cm) :	44.5 x 43 x 34
Per 20 feet container :	530 Carton



Details :

Product Name :	Mini Stitch Cutter with Handle
Size of Inner Box (cm) :	14.5 x 3.5 x 7
Unit/Box :	10 Pcs
Outer Boxes :	--
Export carton :	--
Export Carton size (cm) :	--
Per 20 feet container :	--



Chisel Blades

Fine Stainless Steel

Description

FINE (CHISEL) Surgical Blades Stainless Steel are available in the following sizes:
Nos. T-SM61, T-SM61V, T-SM62, T-SM62V, T-SM63, T-SM64, T-SM65, T-SM66, T-SM67, T-SM68
T-SM69, T-SM69B, T-SM90, T-SM91. Each blade is protected by VCI paper in a peel-open
foil sachet. Sterilized by Gamma Radiation of 2.5 M. Rads.

Details :

Product Name : Chisel Blades (Fine)
Size of Inner Box (cm) : 82 x 32 x 22
Unit/Box : 25 pcs
Outer Boxes : 25 Boxes x 25 Units
Export carton : 25 Boxes x 10 Outer
Export Carton size (cm) : 43.5 x 39 x 38.5
Per 20 feet container : 420 Carton

Size	Product Code	Handle
61	TSM61	TSMH
61V	TSM61V	TSMH
62	TSM62	TSMH
62V	TSM62V	TSMH
63	TSM63	TSMH
64	TSM64	TSMH
65	TSM65	TSMH
66	TSM66	TSMH
67	TSM67	TSMH
68	TSM68	TSMH
69	TSM69	TSMH
69B	TSM69B	TSMH
90	TSM90	TSMH
91	TSM91	TSMH



Gouge Blades & Handles

Carbon Steel

Description

These surgical blades are used for feet surgeries and podiatry purposes i.e. manicure and pedicure. Gouge Blades are used by various podologists for hand, feet and skin care. They are available in the following sizes: Nos. 1, 2, 3, 4, 5, 6, 8, 10, 12, 15 and 18.

Each blade is protected by VCI paper in a peel- open foil sachet. Gouge blades fit into handle no. 3 and pedal gouge handles.



Size	Product Code	Fitment No.	Handle
1	TGCS1	3	TGMH
2	TGCS2	3	TGMH
3	TGCS3	3	TGMH
4	TGCS4	3	TGMH
5	TGCS5	3	TGMH
6	TGCS6	3	TGMH
8	TGCS8	3	TGMH
10	TGCS10	3	TGMH
12	TGCS12	3	TGMH
15	TGCS15	3	TGMH
18	TGCS18	3	TGMH

Details :

Product Name : Gouge Blades No. 1 to 5
Size of Inner Box (cm) : 10 x 10 x 3.3
Unit/Box : 100 pcs
Outer Boxes : 10 Boxes x 100 Units
Export carton : 10 Boxes x 10 Outer
Export Carton size (cm) : 59.5 x 48 x 21
Per 20 feet container : 490 Carton

Details :


Product Name : Gouge Blades No. 6 to 18
Size of Inner Box (cm) : 13.8 x 10 x 3.3
Unit/Box : 100 pcs
Outer Boxes : 10 Boxes x 100 Units
Export carton : 10 Boxes x 10 Outer
Export Carton size (cm) : 57 x 39 x 30
Per 20 feet container : 400 Carton

Disposable Scalpels

Stainless Steel Blade

Features

- Pre-assembled stainless steel surgical blades mounted on ABS plastic handle.
- The device is used to prevent injury during fixing the blade on the metal handle.
- Fluted finger grips and the blades are safely protected by LDPE guard.
- Individually packaged in a peel pouch.
- Convenience of single-use.
- Size clearly marked on each scalpel blade.
- Each dispenser box contains 10 scalpels.
- Shelf life of 5 years



Size	Product Code	Fitment No.	Size	Product Code	Fitment No.
9	TDSS9	3	18	TDSSI8	4
10	TDSSI0	3	19	TDSSI9	4
10A	TDSSI0A	3	20	TDSS20	4
11	TDSSI1	3	21	TDSS21	4
11K	TDSSI1K	3	22	TDSS22	4
12	TDSSI2	3	22A	TDSS22A	4
13	TDSSI3	3	23	TDSS23	4
14	TDSSI4	3	24	TDSS24	4
15	TDSSI5	3	25	TDSS25	4
15B	TDSSI5B	3	26	TDSS26	4
15C	TDSSI5C	3	36	TDSS36	4
16	TDSSI6	3			
17	TDSSI7	3			

Details :

Product Name :	Disposable Scalpels
Size of Inner Box (cm) :	20 x 7 x 3.5
Unit/Box :	10 pcs
Outer Boxes :	10 Boxes x 10 Units
Export carton :	10 Boxes x 10 Outer
Export Carton size (cm) :	43.5 x 39 x 38.5
Per 20 feet container :	420 Carton

Universal Safety Scalpel

Stainless Steel Blade

Features

- **Safety Mode** : When you retract back after fully extending the blade, it makes a click sound to denote it is in safety mode. Safety mode means that the scalpel is not fully extended. This feature avoids injuries when the blade is half way extended as it retracts back and does not pierce through. It has to be fully extended to use.
- **Safety Lock**: if you completely retract back it makes a loud click denoting it is in safety lock position. Safety lock position means that the blade cannot be reused once it is in lock position. This feature forces the doctor to dispose the product and minimizes the risk of transmissibility of blood-borne infectious agents.
- Compatible for a left handed as well as right handed surgeon.
- Safety feature in pre use and post use of product.
- One-handed activation allowing the clinician to focus on the surgical site.
- The retractable design protects the blade from exposure, minimizing the risk of accidental cuts.
- Universal size handle that fits all types of hand sizes, and ideal for kits, trays & packs.
- One handle size fits all blade sizes.



Size	Product Code	Fitment No.	Size	Product Code	Fitment No.
9	TSF9	3	18	TSF18	4
10	TSF10	3	19	TSF19	4
10A	TSF10A	3	20	TSF20	4
11	TSF11	3	21	TSF21	4
11K	TSF11K	3	22	TSF22	4
13	TSF13	3	22A	TSF22A	4
14	TSF14	3	23	TSF23	4
15	TSF15	3	24	TSF24	4
15B	TSF15B	3	25	TSF25	4
15C	TSF15C	3			
16	TSF16	3			
17	TSF17	3			

Details :

Product Name :	Safety Scalpel
Size of Inner Box (cm) :	17 x 90 x 4
Unit/Box :	10 pcs
Outer Boxes :	5 Boxes x 10 Units
Export carton :	5 Boxes x 20 Outer
Export Carton size (cm) :	54 x 39 x 44.5
Per 20 feet container :	280 Carton


Biopsy Punch

Stainless Steel Blade

Description

Biopsy punch is a sharp rounded blade tip used to remove tissue to be tested during procedures. The instrument is sterile with a round plastic handle. The punch is available in a variety of tip sizes and is sold in both sterile single pack as well as in bulk quantities.

Finely ground seamless, stainless steel blade ensures safe, clean and easy sample extraction. The sleek ribbed handle clear view of the biopsy site, comfort and control. Virgin grade ABS plastic provides strength and sturdiness.



Size	Product Code	Product Code
1.0	TBPS1.0	TBPP1.0
1.5	TBPS1.5	TBPP1.5
2.0	TBPS2.0	TBPP2.0
2.5	TBPS2.5	TBPP2.5
3.0	TBPS3.0	TBPP3.0
3.5	TBPS3.5	TBPP3.5
4.0	TBPS4.0	TBPP4.0
5.0	TBPS5.0	TBPP5.0
6.0	TBPS6.0	TBPP6.0
7.0	TBPS7.0	TBPP7.0
8.0	TBPS8.0	TBPP8.0

Details :

Product Name :	Biopsy Punch
Size of Inner Box (cm) :	17.2 x 9 x 3.6
Unit/Box :	10 pcs
Outer Boxes :	10 Boxes x 10 Units
Export carton :	10 Boxes x 5 Outer
Export Carton size (cm) :	40 x 20 x 52.5
Per 20 feet container :	660 Carton

Dermal Curette

Stainless Steel Blade

Description

The curette has a firm place in the range of dermatological, diagnostic and therapeutic instruments. It is one of the most versatile instruments. For that reason, it is increasingly also used in other disciplines, such as plastic surgery and gynaecology. The aim of a curette is the removal of a superficial skin changes with minimal damage to the surrounding tissue. The curette is extremely sharp thus ensuring precise and gentle tissue removal. It can also be used for the treatment of superficial malignant tumours. The outstanding cutting ability makes working with it much easier compared with conventional products.



DCS



DCP

Size	Product Code	Product Code
2.0	TDCS2.0	TDCP2.0
3.0	TDCS3.0	TDCP3.0
4.0	TDCS4.0	TDCP4.0
5.0	TDCS5.0	TDCP5.0
7.0	TDCS7.0	TDCP7.0

Details :

Product Name :	Dermal Curette
Size of Inner Box (cm) :	21 x 7 x 4
Unit/Box :	10 pcs
Outer Boxes :	5 Boxes x 10 Units
Export carton :	5 Boxes x 5 Outer
Export Carton size (cm) :	46.5 x 46 x 25
Per 20 feet container :	500 Carton

Ophthalmic Knives

Slit-Blades (Keratome) For Phaco Stab Incision



Product No. Description

P-912301 2.3 MM Sharp Tip Straight
P-912501 2.5 MM Sharp Tip Straiaht
P-912601 2.65 MM Sharp Tip Straight
P-912801 2.8 MM Sharp Tip Straiaht
P-912901 3.0 MM Sharp Tip Straight
P-913201 3.2 MM Sharp Tip Straiaht
P-913501 3.5 MM Sharp Tip Straight
P-915101 5.1 MM Sharp Tip Straiaht

Product No. Description

P-912361 2.3 MM Sharp Tip Angled 45° 19 Gauge
P-912561 2.5 MM Sharp Tip Angled 45° 19 Gauge
P-912661 2.65 MM Sharp Tip Angled 45° 19 Gauge
P-912861 2.8 MM Sharp Tip Angled 45° 19 Gauge
P-912961 3.0 MM Sharp Tip Angled 45° 19 Gauge
P-913261 3.2 MM Sharp Tip Angled 45° 19 Gauge
P-913561 3.5 MM Sharp Tip Angled 45° 19 Gauge
P-915161 5.1 MM Sharp Tip Angled 45° 19 Gauge

Double Bevel

P-912808 2.8 MM Sharp Tip Straight
P-912908 3.0 MM Sharp Tip Straight
P-913208 3.2 MM Sharp Tip Straight

Double Bevel

P-912868 2.8 MM Sharp Tip Angled 45° 19 Gauge
P-912968 3.0 MM Sharp Tip Angled 45° 19 Gauge
P-913268 3.2 MM Sharp Tip Angled 45° 19 Gauge

Lance Tip Blades For Initial Incision



Product No. Description

P-931501 15° Degree 24 Gauge
P-933001 30° Degree 24 Gauge
P-934501 45° Degree 24 Gauge

MVR Blades For Water Tight Self Sealing Incision



Product No. Description

P-975559 MVR 24 Gauge
P-975560 MVR 20 Gauge
P-975561 MVR 19 Gauge

Product No. Description

AC STILETIO BLADES
For Blumenthal ar Side Part Incision
P-985560 MVR 20 Gauge Angled 45°
P-985561 MVR 19 Gauge Angled 45°

Slit-Blades (Keratome) For I.O.L. Enlargers



Product No. Description

P-914001 4.0 MM Sharp Tip Straight
P-914101 4.1 MM Sharp Tip Straight
P-915101 5.1 MM Sharp Tip Straight
P-915201 5.2 MM Sharp Tip Straight
P-915501 5.5 MM Sharp Tip Straight
P-916001 6.0 MM Sharp Tip Straight
P-916201 6.2 MM Sharp Tip Straight

Product No. Description

P-914061 4.0 MM Sharp Tip Angled 45° 16 Gauge
P-914161 4.1 MM Sharp Tip Angled 45° 16 Gauge
P-915161 5.1 MM Sharp Tip Angled 45° 16 Gauge
P-915261 5.2 MM Sharp Tip Angled 45° 16 Gauge
P-915561 5.5 MM Sharp Tip Angled 45° 15 Gauge
P-916061 6.0 MM Sharp Tip Angled 45° 15 Gauge
P-916261 6.2 MM Sharp Tip Angled 45° 15 Gauge

Crescent Blades For Groove Or Pocket Or Tunnel Incision



Product No. Description

P-950001 Crescent Sharp Tip Straight

Double Bevel

P-950001 Crescent Sharp Tip Straight

Product No. Description

P-950002 Crescent Bevel Up Angled 45° 20 Gauge
P-950003 Crescent Bevel Down Angled 45° 20 Gauge
P-950004 Crescent Angled 45° 20 Gauge

Details :

Product Name : Ophthalmic Knives
Size of Inner Box (cm) : 20 x 5 x 4
Unit/Box : 6 pcs
Outer Boxes : 25 Boxes x 6 Units
Export carton : 25 Boxes x 2 Outer
Export Carton size (cm) : 44 x 28.5 x 23
Per 20 feet container : 1000 Carton



1, L.S.C., Okhla Industrial Area, Phase-II, New Delhi-110 020, India
Tel +91-11-46436600, 46436601, 46436666, 26389812, 13
Fax +91-11-26389815, 41616555
sales@paramountblades.com
export2@paramountblades.com
www.paramountsurgimed.com

(Disclaimer- All measurements in weight & size are indicative for reference only.)