




The Paramount Philosophy:

With a vision of being a global leader in the medical disposable industry, Paramount Surgimed believes in providing its customers a timely supply of quality products at a competitive price. The company came into existence in 1993 under the Grover Group of companies capitalizing on the rich and diverse experience of over 66 years. Now, the brand is symbolized with trust and quality in the medical industry. The company started with the manufacturing of Surgical Blades and is now a leading supplier of a huge range of disposable products.

Paramount exports to over 100 countries, has 40 distributors and 400 Dealers across India. Our State-of-the-Art production facility compares to any international standard and is accredited with ISO: 13485:2016, CE MDR Certificate and US FDA Device Listing.

The Paramount Vision:

Paramount Surgimed and all its employees share the vision of their dynamic Managing Director, Mr. Shaily Grover to be a global leader in the Medical Industry. Mr. Grover has ensured that the Group's reputation amongst its associates is synonymous with top quality, updated technology, dedicated approach and customer friendliness.



Corporate Office – Okhla, New Delhi



Factory – Bhiwadi, Rajasthan



Factory – Bawal, Haryana



Factory – Bhiwadi, Rajasthan

Surgical Blades/Scalpel Blade Remover	02 - 03
Disposable Scalpels	04
Safety Scalpel	05 - 07
Gouge Blades & Handles	08
Chisel/Fine/Micro Surgical Blade	09
Stitch Cutters & Thumb Scalpel	10-11
Skin Graft Blade	12
Myringotomy Knives	13
Biopsy Punch	14 - 15
Dermal Curette	16 - 17
Ophthalmic Knives	18 - 19
Soft Core Biopsy (Automatic)	20
Bone marrow Biopsy Needle (Manual)	21
Hernia Mesh	22
Skin Stapler/Skin Staple Remover	23
Sutures	24 - 25
Foley Catheter	26 - 27
Microtome Blade	28
Blood Lancet, Twisted & Flyfix	29
Disposable Adult Diapers	30 - 31
Under Pad - AroKleen & PupPee	32
Disposable Prep Razor	33
Face Mask	34 - 35
Wipes	36 - 38
Lung Exerciser, Hot Water Bottle	39
Infrared Thermometer, Pulse Oximeter & Air Mattress	40
Certificates	41

Surgical Blades

Carbon Steel & Stainless Steel

Features

- Stainless Steel Surgical Blades and Carbon Steel Surgical Blades.
- Packed in peel open pouches with a VCI protective inner paper.
- Shelf life of 5 years - Sterile Product
- Stainless Steel (Grade F as per BSEN ISO 7153-1:2001)
- Carbon Steel (Grade- As per BS 2982:1992)
- Complies the Hardness requirement of BS 2982:1992 standard.
- Complies the fitment requirement of ISO 7740:1985 standard requirement.



FITMENT 3



FITMENT 4



Details:

Product Name:	Surgical Blades
Size of Inner Box (cm):	8.5 x 7.0 x 2.7
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	42.5 x 39 x 17
Per 20 feet container:	825 Carton

Post Mortem Blades –Size 40

Features:

- The PM40 range of blades is designed for use with screw type heavy duty handles.
- Heavy duty blade.
- Box of 10.



Post Mortem Blades –Size 60

Features:

- Used in the fields of Pathology and Histopathology for autopsy.
- It is used with handle No.8.
- Available in both material carbon steel and stainless steel.
- Box of 50.
- Individually packed in easy to open peel pouch.



Scalpel Blade Remover

Description

The Blade Remover enables safe, fast and easy removal of used surgical blades from any size or type of blade handles. It improves safety and hygiene by eliminating direct contact with blades and facilitates replacement.



Size	Product Code (Carbon Steel)	Product Code (Stainless Steel)	Fitment No
9	PCS9	PSS9	3
10	PCS10	PSS10	3
10A	PCS10A	PSS10A	3
11	PCS11	PSS11	3
11K	PCS11K	PSS11K	3
11P	PCS11P	PSS11P	3
12	PCS12	PSS12	3
12D	PCS12D	PSS12D	3
13	PCS13	PSS13	3
14	PCS 14	PSS14	3
15	PCS15	PSS15	3
15B	PCS15B	PSS15B	3
15C	PCS15C	PSS15C	3
15D	PCS15D	PSS15D	3
15S	PCS15S	PSS15S	3
16	PCS16	PSS16	3
17	PCS17	PSS17	3
18	PCS18	PSS18	4
19	PCS19	PSS19	4
20	PCS20	PSS20	4
21	PCS21	PSS21	4
22	PCS22	PSS22	4
22A	PCS22A	PSS22A	4
23	PCS23	PSS23	4
24	PCS24	PSS24	4
24D	PCS24D	PSS24D	4
25	PCS25	PSS25	4
26	PCS26	PSS26	4
36	PCS36	PSS36	4
36D	PCS36D	PSS36D	4
40	PCS40	PSS40	–
40B	PCS40B	PSS40B	-
60	PCS60	PSS60	8
60B	PCS60B	PSS60B	8


Product Name:	Scalpel Blade Remover
Material:	Plastic Box
Capacity:	100 Pieces of Blades
Inner Box(cm):	17X12.5X5.8
Unit:	10 Pcs
Export Carton:	10 Boxes
Export Size(cm):	36X31X14.5

Disposable Scalpels

Stainless Steel blade

Features

- Pre-assembled stainless steel surgical blades mounted on ABS plastic handle.
- The device is used to prevent injury during fixing the blade on the metal handle.
- Fluted finger grips and the blades are safely protected by LDPE Guard / Cap.
- Individually packed in easy to open peel pouch.
- Single use design for maximum safety and convenience.
- Size clearly marked on each scalpel blade and blister pouch.
- Each dispenser box contains 10 scalpels.
- Shelf life of 5 years - Sterile Product.



Details:

Product Name:	Disposable Scalpel
Size of Inner Box (cm):	21 x 7 x 3
Unit/Box:	10 Pcs
Export carton:	100 Boxes
Export Carton size (cm):	44 x 37 x 33
Per 20 feet container:	370 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
9	PDSS9	3	18	PDSS18	4
10	PDSS10	3	19	PDSS19	4
10A	PDSS10A	3	20	PDSS20	4
11	PDSS11	3	21	PDSS21	4
11K	PDSS11K	3	22	PDSS22	4
12	PDSS12	3	22A	PDSS22A	4
13	PDSS13	3	23	PDSS23	4
14	PDSS14	3	24	PDSS24	4
15	PDSS15	3	24D	PDSS24D	4
15D	PDSS15D	3	25	PDSS25	4
15S	PDSS15S	3	36	PDSS36	4
15B	PDSS15B	3	11P	PDSS11P	3
15C	PDSS15C	3	12D	PDSS12D	4
16	PDSS16	3	34	PDSS34	4
17	PDSS17	3	36D	PDSS36D	4



Safety Scalpel

Stainless Steel blade

Features

- It is a safe alternative to the traditional scalpel and fulfills the requirement to reduce sharp injuries.
- Ribbed on both side of the sheath which provide extra gripping for left handed as well as right handed surgeon.
- One-handed activation allowing the surgeon to focus on the surgical site.
- Audible clicks makes the scalpel more secure and protective.
- "Used & Un-Used" indications prevent contamination and transmission of blood borne disease.
- Permanent safety lock prevents accidental cuts.
- Latex Free.
- Measurement scale 0 to 5 cm is embossed on the handle.
- Ergonomically designed safety scalpel.
- ABS plastic handle and sheath gives extra strength to the product.
- Universal size handle that is ideal for kits, trays & packs.

Details:

Product Name:	Safety Scalpel
Size of Inner Box (cm):	21 x 7 x 3
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	44 x 37 x 17
Per 20 feet container:	300 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
9	PNSF9	3	16	PNSF16	4
10	PNSF10	3	17	PNSF17	4
10A	PNSF10A	3	18	PNSF18	4
11	PNSF11	3	19	PNSF19	4
11K	PNSF11K	3	20	PNSF20	4
11P	PNSF11P	3	21	PNSF21	4
12	PNSF12	3	22	PNSF22	4
12D	PNSF12D	3	22A	PNSF22A	4
13	PNSF13	3	23	PNSF23	4
14	PNSF14	3	24	PNSF24	4
15	PNSF15	3	25	PNSF25	4
15B	PNSF15B	3	36	PNSF36	4
15C	PNSF15C	3	24D	PNSF24D	4
15D	PNSF15D	3	34	PNSF34	4
15S	PNSF15S	3	36D	PNSF36D	4



Safety Scalpel

Stainless Steel blade

Features

- Safe alternative to traditional scalpels, reducing the risk of sharp injuries.
- Dual-sided ribbed sheath provides superior grip for both left- and right-handed surgeons.
- One-handed activation allowing the surgeon to maintain focus on the surgical site.
- Audible click ensures secure and protective blade engagement.
- Equipped with a permanent safety lock for easy disposal and enhanced protection against accidental cuts.
- Latex-free construction.
- Embossed measurement scale (0–5 cm) on the handle for precision.
- Ergonomically designed for comfort & control with smooth sliding lock.
- Durable ABS plastic handle and sheath provide extra strength to the product.
- Standard handle size-ideal for surgical kits, trays, and procedure packs.

Details:

Product Name:	Safety Scalpel
Size of Inner Box (cm):	21 x 7 x 3.5
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	44 x 37 x 17
Per 20 feet container:	300 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
10	PNSS10	3	20	PNSS20	4
11	PNSS11	3	21	PNSS21	4
12	PNSS12	3	22	PNSS22	4
13	PNSS13	3	23	PNSS23	4
15	PNSS15	3	24	PNSS24	4
15C	PNSS15C	3			



Universal Safety Scalpel

Stainless Steel blade

Features

- **Safety Mode:** When you retract back after fully extending the blade, it makes a click sound to denote it is in safety mode. Safety mode means that the scalpel is not fully extended. This feature avoids injuries when the blade is half way extended as it retracts back and does not pierce through. It has to be fully extended to use.
- **Safety Lock:** if you completely retract back it makes a loud click denoting it is in safety lock position. Safety lock position means that the blade cannot be reused once it is in lock position. This feature forces the doctor to dispose the product and minimizes the risk of transmissibility of blood-borne infectious agents.
- Compatible for a left handed as well as right handed surgeon.
- Safety feature in pre use and post use of product.
- One-handed activation allowing the clinician to focus on the surgical site.
- The retractable design protects the blade from exposure, minimizing the risk of accidental cuts.
- Universal size handle that fits all types of hand sizes, and ideal for kits, trays & packs.
- One handle size fits all blade sizes.

Details:

Product Name:	Universal Safety Scalpel
Size of Inner Box (cm):	17 x 8 x 3.5
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	43 x 37 x 20.5
Per 20 feet container:	280 Carton

Size	Product Code	Fitment No	Size	Product Code	Fitment No
9	PSF9	3	11P	PSF11P	3
10	PSF10	3	18	PSF18	4
10A	PSF10A	3	19	PSF19	4
11	PSF11	3	20	PSF20	4
11K	PSF11K	3	21	PSF21	4
12	PSF12	3	22	PSF22	4
13	PSF13	3	22A	PSF22A	4
14	PSF14	3	23	PSF23	4
15	PSF15	3	24	PSF24	4
15B	PSF15B	3	25	PSF25	4
15C	PSF15C	3	36	PSF36	4
15D	PSF15D	3	12D	PSF12D	4
15S	PSF15S	3	24D	PSF24D	4
16	PSF16	3	36D	PSF36D	4
17	PSF17	3	34	PSF34	4




Gouge Blades

Carbon Steel & Stainless Steel

Description

These surgical blades are used for foot surgeries and podiatry purposes. Gouge Blades are used by various podiatrist for hand, feet and skin surgery. They are available in the following sizes: Nos. 1, 2, 3, 4, 5, 6, 8, 10, 12, 15 and 18. Each blade is protected by VCI paper in a peel- open foil sachet. Gouge blades fit into handle no. 3 and special gouge handles.



Details:

Product Name:	Gouge Blades No. 1 to 5
Size of Inner Box (cm):	10 x 10 x 3.3
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	59.5 x 47 x 19.5
Per 20 feet container:	425 Carton

Size	Product Code (CS)	Product Code (SS)	Size	Product Code	Product Code	Fitment No	Handle
1	PGCS1	PGSS1	1V	PGCS1V	PGSS1V	3	PGMH
2	PGCS2	PGSS2	2V	PGCS2V	PGSS2V	3	PGMH
3	PGCS3	PGSS3	3V	PGCS3V	PGSS3V	3	PGMH
4	PGCS4	PGSS4	1R	PGCS1R	PGSS1R	3	PGMH
5	PGCS5	PGSS5	2R	PGCS2R	PGSS2R	3	PGMH
6	PGCS6	PGSS6	3R	PGCS3R	PGSS3R	3	PGMH
8	PGCS8	PGSS8				3	PGMH
10	PGCS10	PGSS10				3	PGMH
12	PGCS12	PGSS12				3	PGMH
15	PGCS15	PGSS15				3	PGMH

Details:

Product Name:	Gouge Blades No. 6 to 18
Size of Inner Box (cm):	13.8 x 10 x 3.3
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	57 x 39 x 30
Per 20 feet container:	390 Carton



Chisel Micro Surgical Blades

Fine Stainless Steel

Description

Chisel Blade is made of Stainless steel (Grade F) and Carbon steel (As per BS 2982). The Blade is packed in peel apart Aluminum Pouch with VCI liner to protect the blade from corrosion. The device is sterilized by Gamma Radiation & Expiry 5 years from the date of manufacturing. The device is used for incision/ cuts during surgery for all blade sizes. The device is under the category of surgically invasive device and for transient use. The device is contact with intact skin and mucosal membrane.

Details:

Product Name:	Chisel Micro Surgical Blades
Size of Inner Box (cm):	8.2 x 3.2 x 2.2
Unit/Box:	25 Pcs
Outer Boxes:	25 Boxes x 25 Units
Export carton:	25 Boxes x 10 Outer
Export Carton size (cm):	47.5 x 39 x 17
Per 20 feet container:	825 Carton

Size	Product Code	Size	Product Code	Handle
61	PSM61	61S	PSM61S	PSMH
62	PSM62	62S	PSM62S	PSMH
63	PSM63	65A	PSM65A	PSMH
64	PSM64			PSMH
65	PSM65			PSMH
67	PSM67			PSMH
68	PSM68			PSMH
69	PSM69			PSMH
90	PSM90			PSMH
91	PSM91			PSMH



Stitch Cutters

Stainless Steel & Carbon Steel

Description

Stitch Cutter is used for general suture removal. Appearing similar to a scalpel, this instrument requires no handle and is a simple and convenient, way to remove basic sutures. The slightly curved, flat, bull nose tip slides under the suture and when rotated, gently removes the suture from the site. Long stitch cutter is made of stainless steel and short stitch cutter is made of carbon steel. The mini stitch cutter is made of stainless steel can be used with handle No.3L, 7 & 9. Alternatively, it is also available on pre-mounted ABS disposable plastic handle.

Features

- 5 Year Shelf Life – Sterile Product.
- Each blade is protected by VCI paper in a peel-open foil sachet.
- Convenience of single-use disposables.
- Sizes clearly marked on each blade.
- Each box contains 100 pcs.

Size	Product Code (Carbon Steel)	Product Code (Stainless Steel)	Handle No.
LSC	–	PLSCSS	-
SSC	PSSCCS	–	-
MSC	–	PMSCSS	PCS3L

Details:

Product Name:	Mini Stitch Cutter
Size of Inner Box (cm):	8.5 x 7.0 x 2.7
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	42.5 x 39 x 17
Per 20 feet container:	825 Carton



Details:

Product Name:	Short Stitch Cutter
Size of Inner Box (cm.):	11 x 8 x 3.5
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm.):	46.5 x 40.0 x 25.5
Per 20 feet container:	550 Carton



Long Stitch Cutters

Stainless Steel & Carbon Steel

Details:

Product Name:	Long Stitch Cutter
Size of Inner Box (cm):	14.8 x 7 x 3.1
Unit/Box:	100 Pcs
Outer Boxes:	10 Boxes x 100 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	45 x 34 x 36
Per 20 feet container:	475 Carton



Mini Stitch Cutter with Handle

Stainless Steel

Details:

Product Name:	Mini Stitch Cutter with Handle
Size of Inner Box (cm):	14.5 x 70 x 3.5
Unit/Box:	Bulk Pack of 50 Pcs
Outer Boxes:	Not Using
Export carton:	10 Bags of 50 Pcs x 10 Outer
Export Carton size (cm):	42 x 32 x 39
Per 20 feet container:	290 Carton



Thumb Scalpel

Stainless Steel


Description

Disposable thumb scalpel or pocket scalpel is used in place of a surgical blade. This reduces or even prevents the incidences of injuries that can occur during the process of fitting or removing contaminated blades from handles.

The disposable thumb scalpels have a protective cap/guard for the blade that can be used so the scalpel can be safely disposed in the operation theater/room disposal container.

Details:

Product Name:	Thumb Scalpel
Product Code:	PTSS
Size of Inner Box (cm):	15.5x 11.5 x 8
Unit/Box:	25 Pcs
Export carton:	20 Boxes x 25 Units
Export Carton size (cm):	60 x 33.5 x 21
Per 20 feet container:	700 Carton



Skin Graft Blade

Carbon Steel

Description


Skin Graft Blade - Simplex, Duplex

Surgical Skin Graft Blades are 158 mm long, individually wrapped in VCI paper in single peel packets in packs of 10. These blades can be used with Watson, Braithwaite or Cobbett Knives.

Details:

Product Name:	Skin Graft Blades
Size of Inner Box (cm):	20.7 x 3.6 x 2
Unit/Box:	10 Pcs
Outer Boxes:	10 Boxes x 10 Units
Export carton:	10 Boxes x 10 Outer
Export Carton size (cm):	44 x 24.5 x 22.5
Per 20 feet container:	1100 Carton

Size	Carbon Steel	Stainless Steel
SG	PCSG	PSSG
SGS	PCSGS	PSSGS
SGD	PCSGD	PSSGD



Metal Handle

Description

Metal handles comes in different sizes and are used according to the applications. The most popular stainless steel handle sizes used are size #3 and size #4 handles.

Gouge Handle - PGMH



Chisel Handle - PSMH



Metal Handle Fit.3 - PSSH3



Metal Handle Fit.4 - PSSH4



Myringotomy Knives

Description

Myringotomy is a surgical procedure of the eardrum or tympanic membrane. The procedure is performed by making a small incision with a myringotomy knife through the layers of tympanic membrane.

Features:

- Optimised micrometre-accurate precision grinding.
- Double edged allowing for cutting movements in both directions.
- Ultra-sharp cutting edges of superb durability.
- Consistently sharp, angled blade facilitates a clean cut.
- Spear tip, Narrow shaft.
- Lance tip, Downward cut, Narrow shaft .
- With and Without Handle.



Details:

Product Name:	Myringotomy Knives
Size of Inner Box (cm):	24.2 x 8 x 2.5
Unit/Box:	5 Pcs
Outer Boxes:	--
Export carton:	50 Boxes
Export Carton size (cm):	50.5 x 42 x 14
Per 20 feet container:	1000 Carton

Size	Product Code
Lance	PSSML
Spear	PSSMS



Suction Needle

Suction Tube

Description

The suction controller tips are a range of instruments used in aural micro suction. It connects to an micro suction machine through a fenestrated suction controller. Each one is colour coded for easy identification of sizes and features a sterilized a stainless steel tube. Each suction controller tip is individually sealed. It is a disposable instrument intended for single use

Features

- Angled Cannula.
- Luer Conical connection Plug.
- Rounded Cannula Point.
- Cannulas can be use for dissection of soft issue.



Details:

Product Name:	Suction Needle
Size of Inner Box small (cm):	18x7x3
Size of Inner Box large (cm):	22x7x3
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size small (cm):	38x37x17.5
Export Carton size large (cm):	46x37x17.5

Metal Suction Needle
2.10 mm x 70mm

Colour	Green
Size	14G-2.10MM X 70MM

Metal Suction Needle
60mm Length

Colour	Size
White	17G-1.47MM X 60mm
Brown	19G-1.07MM X 60mm
Black	22G-0.70MM X 60mm
Purple	27G-0.40MM X 60mm

Micro Suction Tube

Colour	Size
Blue	3.00mm X 130mm
Yellow	4.00mm X 130mm

Micro Suction Tube
80mm Length

Colour	Size
Brown	1.00mm X 80mm
White	1.50mm X 80mm
Green	2.00mm X 80mm
Orange	2.50mm X 80mm
Yellow	4.00mm X 80mm

Biopsy Punch

Stainless Steel blade

Description

Biopsy punch is a sharp rounded blade tip used to remove tissue to be tested during procedures. The instrument is sterile with a round plastic handle. The punch is available in a variety of tip sizes and is sold in both sterile single pack as well as in bulk quantities.


Finely ground seamless, stainless steel blade ensures safe, clean and easy sample extraction. The sleek ribbed handle clear view of the biopsy site, comfort and control.

Virgin grade ABS plastic provides strength and sturdiness.

Details:

Product Name:	Biopsy Punch
Size of Inner Box (cm):	17 x 8 x 3.5
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	43 x 37 x 20.5
Per 20 feet container:	635 Carton

Size mm	Product Code	Product Code
1.0	BPS1.0	BPP1.0
1.5	BPS1.5	BPP1.5
2.0	BPS2.0	BPP2.0
2.5	BPS2.5	BPP2.5
3.0	BPS3.0	BPP3.0
3.5	BPS3.5	BPP3.5
4.0	BPS4.0	BPP4.0
5.0	BPS5.0	BPP5.0
6.0	BPS6.0	BPP6.0
7.0	BPS7.0	BPP7.0
8.0	BPS8.0	BPP8.0



BPS Features

- Ergonomically designed for precise skin excisions.
- Unique, ribbed, easy-grip barrel and pressure assistance platform.
- Sterile and ready to use.
- Wide range of sizes.



BPP Features

- Premium quality punch.
- Hard blister packaging.
- Option to print custom brand on handle.
- Pressure assistance platform and ribbed handle.
- Micro grinding finish.



Dermal Blade


Stainless Steel blade

Description

The device is intended for cutaneous surgery giving smooth clean excisions with less tissue trauma. It is used to remove surface protuberances, flat lesions, or deep lesions by saucer shaped incision. Dermal blade is made of stainless-steel blade mounted on flexible plastic holder. The device is under category of surgically invasive device and for transient use. It is covering under surface device (nature of body contact) and contact with intact skin and mucosal membrane. The device is invasive device and for transient use.

Details:

Product Name:	Dermal Blade
Product Code:	PDB
Size of Inner Box (cm):	15.5 x 11.5 x 11
Unit/Box:	20 Pcs
Export carton:	20 Boxes
Export Carton size (cm):	60 x 33.5 x 24
Per 20 feet container:	700 Carton



Dermal Curette

Stainless Steel blade

Description

The curette has a firm place in the range of dermatological, diagnostic and therapeutic instruments. It is one of the most versatile instruments. For that reason, it is increasingly also used in other disciplines, such as plastic surgery.


The aim of a curette is the removal of a superficial skin changes with minimal damage to the surrounding tissue.

The curette is extremely sharp thus ensuring precise and gentle tissue removal. The outstanding cutting ability makes working with it much easier compared with conventional products.

Details:

Product Name:	Dermal Curette
Size of Inner Box (cm):	21 x 7 x 3
Unit/Box:	10 Pcs
Export carton:	50 Boxes
Export Carton size (cm):	44 x 37 x 17
Per 20 feet container:	475 Carton

Size mm	Product Code	Product Code	Product Code	Product Code
2.0	DCS2.0	DCP2.0	DCK2.0	DDC2.0
3.0	DCS3.0	DCP3.0	DCK3.0	DDC3.0
4.0	DCS4.0	DCP4.0	DCK4.0	DDC4.0
5.0	DCS5.0	DCP5.0	DCK5.0	DDC5.0
7.0	DCS7.0	DCP7.0	DCK7.0	DDC7.0



DCS

DCS Features


- Seamless stainless steel sharp cutting edge.
- Easy dynamic gripping handle.
- Sterile and ready to use.
- Sharp side indicator.
- Wide range of sizes.



DCK

DCK Features

- Sturdy & Durable Stainless Steel.
- Symmetrically Ground Beveled Edge.
- Ultra Sharp Inside Grinding.
- Sharp side indicator.
- Ribbed on all sides.
- Triangular Handle.
- Sterile and ready to use.
- Wide range of sizes.



DCP

DCP Features


- Long 'Pen Type' Grip.
- Branding Labeling area.
- Ribbed handle for easy control.
- Minimum damage to surrounding tissues.
- Sharp side indicator.
- Sterile and ready to use.
- Wide range of sizes.



DDC

DDC Features

- Seamless stainless steel sharp cutting edge.
- One in metal with a sharp edge.
- Sharp side indicator.
- Covered by a protective cap.
- Ultra Sharp Inside Grinding.
- Sterile and ready to use.
- Wide range of sizes.



Ophthalmic Knives

Description

High-precision surgical knives designed for ophthalmic procedures, ensuring smooth, accurate, and consistent incisions. Available in a wide range of blade geometries and sizes for various surgical needs.

- **Slit / Keratome Knives:** For phaco and IOL incision procedures; available in blunt or sharp tips, straight or angled designs.
- **Crescent Knives:** Ideal for creating precise grooves, pockets, or tunnel incisions.
- **Lance Tip Knives:** Used for initial corneal entry with sharp, angled tips for controlled penetration.
- **MVR Knives:** Designed for self-sealing and watertight incisions with round-wire construction.
- **Spoon & Scleral Knives:** Specialized for limbal, scleral, and tunnel incisions with narrow or wide tips.

All knives are crafted for superior sharpness, durability, and surgeon comfort — supporting the quest for better healthcare.

Round-Wire Knives



Slit-Blade (Keratome) For Phaco Stab Incision

Product No. Description

P-912301 2.3 MM Sharp Tip Straight
P-912501 2.5 MM Sharp Tip Straight
P-912601 2.65 MM Sharp Tip Straight
P-912801 2.8 MM Sharp Tip Straight
P-912901 3.0 MM Sharp Tip Straight
P-913201 3.2 MM Sharp Tip Straight
P-913501 3.5 MM Sharp Tip Straight
P-915101 5.1 MM Sharp Tip Straight

Double Bevel

P-912808 2.8 MM Sharp Tip Straight
P-912908 3.0 MM Sharp Tip Straight
P-913208 3.2 MM Sharp Tip Straight

Product No. Description

P-912361 2.3 MM Sharp Tip Angled 45° 19 Gauge
P-912561 2.5 MM Sharp Tip Angled 45° 19 Gauge
P-912661 2.65 MM Sharp Tip Angled 45° 19 Gauge
P-912861 2.8 MM Sharp Tip Angled 45° 19 Gauge
P-912961 3.0 MM Sharp Tip Angled 45° 19 Gauge
P-913261 3.2 MM Sharp Tip Angled 45° 19 Gauge
P-913561 3.5 MM Sharp Tip Angled 45° 19 Gauge
P-915061 5.1 MM Sharp Tip Angled 45° 19 Gauge

Double Bevel

P-912868 2.8 MM Sharp Tip Angled 45° 19 Gauge
P-912968 3.0 MM Sharp Tip Angled 45° 19 Gauge
P-913268 3.2 MM Sharp Tip Angled 45° 19 Gauge



Slit-Blade (N-Larger) For I.O.L. Enlargers

Product No. Description

P-914001 4.0 MM Blunt Tip Straight
P-915201 5.2 MM Blunt Tip Straight
P-915501 5.5 MM Blunt Tip Straight
P-916001 6.0 MM Blunt Tip Straight
P-916201 6.2 MM Blunt Tip Straight

Product No. Description

P-914061 4.0 MM Blunt Tip Angled 45° 16 Gauge
P-915561 5.5 MM Blunt Tip Angled 45° 15 Gauge
P-916061 6.0 MM Blunt Tip Angled 45° 15 Gauge
P-916261 6.2 MM Blunt Tip Angled 45° 15 Gauge



Crescent Blade For Groove Or Pocket Or Tunnel Incision

Product No. Description

P-950001 Crescent Sharp Tip Straight

Double Bevel

P-950004 Crescent Angled 45° 20 Gauge
P-950005 Sharp Tip Straight

Product No. Description

P-950002 Crescent Bevel Up Angled 45° 20 Gauge
P-950003 Crescent Bevel Down Angled 45° 20 Gauge



Lance Tip Blade For Initial Incision

Product No. Description

P-931501 15° Degree 24 Gauge

P-933001 30° Degree 24 Gauge

P-934501 45° Degree 24 Gauge

Round-Wire Knives



MVR Blade For Water Tight Self Sealing Incision

Product No. Description

P-975559 MVR 24 Gauge

P-975560 MVR 20 Gauge

P-975561 MVR 19 Gauge

Product No. Description

AC STILETIO BLADES

For Blumenthal or Side Port Incision

P-985560 MVR 20 Gauge Angled 45°

P-985561 MVR 19 Gauge Angled 45°



Miniature Blade Alternative To General Surgery For Precise Incision

Product No. Description

SPOON KNIVES

For Groove or Pocket or Tunnel Incision

Product No. Description

P-6821 3.0 MM Angled 60°

P-6821E 3.0 MM Extended Neck Angled 60°



Miniature Blade Alternative To General Surgery For Precise Incision

Product No. Description

SCLERAL KNIVES

For Limbal and Scleral Incision

P-5700 Wide Tip Blade

P-5710 Narrow Tip Blade



Details:

Product Name:	Ophthalmic Knives
Size of Inner Box (cm):	20 x 5 x 4
Unit/Box:	6 Pcs
Export carton:	25 Boxes x 2 Outer
Export Carton size (cm):	42.5 x 27.5 x 22.5
Per 20 feet container:	1000 Carton


Soft Core Biopsy (Automatic) Biopsy Needle

Features:

- Single handed use allows for faster, easier use with more control with gloved hand.
- Dual triggers allows for flexibility: physician can use different biopsy techniques.
- Compatibility with TruGuide Coaxial Biopsy Needle enhances efficiency and accuracy.
- One-handed cocking and lightweight ergonomic design improve both handling and control.
- Two firing buttons accommodate your preference.
- Smooth surface: Less friction, Consistent performance, frictionless insertion reduce vessel trauma and pain.
- Thin Wall: Better Inner diameter helps to get better Inner diameter of cannula provides more space.

Details:

Product Name:	Soft Core Biopsy (Automatic)
Size of Inner Box (cm):	49X16.5X7
Unit/Box:	5 Pcs
Export carton:	10 Boxes
Export Carton size (cm):	51.5X35X37




Bone Marrow Biopsy (Manual) Biopsy Needle

Features:

- "OPEN" and "CLOSED" locking lever captures and secures intact specimen.
- Extra sharp tip allows for easy penetration of the cortex.
- "Quintuple peaks" on cannula tip to provide superior coring ability.
- Ergonomic handle designed to minimize hand and arm fatigue.
- Marked probe to verify specimen length before extraction.
- Tip probe guide ensures end user safety from accidental needle stick.
- Needle protection included for securing sharps after use.

Details:

Product Name:	Bone Marrow Biopsy (Manual)
Size of Inner Box (cm):	41.5X11.5X6.5
Unit/Box:	5 Pcs
Export carton:	10 Boxes
Export Carton size (cm):	44X25X35



Hernia Mesh

ParaLENE Mesh

Description

- Laser cut providing precise shape and preventing edges from fraying.
- Permeable giving it high capacity for colonization, good tissue ingrowth and integration.
- Transparent and flexible for easy use in coelioscopy and laparotomy.
- Fits perfectly against the cavity wall and avoids fragmentation.
- Excellent multi-directional mechanical properties and exceptional resistance to tearing.
- No shrinkage providing long term material stability.
- High strength and burst resistance for permanent support.
- Good shape memory.


Features

- Abdominal wall reinforcement.
- Treatment of incisional, inguinal, femoral, umbilical and ventral hernias.
- Use in both coelioscopy and laparotomy.



Details :

Product Name	Hernia Mesh
Composition	Polypropylene
Construction	Monofilament
Absorption Time	Non-Absorbable
Sterilization	100% Ethylene Oxide
Shelf Life	5 years
Inner Box(cm)	27X24X0.8
Unit	6 Pcs
Export Carton	20 Boxes
Export Carton Size(cm)	29x26x18


	Paralene Mesh	Paralene Lite	Paralene Ultra Lite
Mesh Weight	80 g/m ²	45 g/m ²	28 g/m ²
Mesh Thickness	0.50 mm	0.43 mm	0.30 mm
Pore Size	0.5 x 0.7 mm	1.3 x 1.0 mm	1.0 x 1.0 mm
Elongation	89%	77%	51%
Maximum Force	25 kgf	19 kgf	15 kgf



Paralene Mesh
(magnified view)



Paralene Lite
(magnified view)



Paralene Ultra Lite
(magnified view)

Skin Stapler

ParaSTAPLER

Description

The Skin Stapler is a sterile, single patient use instrument designed to deliver rectangular, stainless steel staples for routine wound closure. The staples are made of lubricant-coated stainless steel. The Skin Stapler places single staples to close surgical incisions. It is supplied sterile and is preloaded for single patient use.

Features

- Easy staple extraction.
- Superior cosmetic results - promotes blood supply and faster scar-free healing.
- Minimizes staple rotation.
- Uncompromised visibility, comfort and control.
- Staples are coated with lubricant.
- 35 Clips.
- Easy Operation.
- Shorten suturing time.

Details:

Product Name:	Skin Stapler
Staple Dimensions:	Wide
Wire Diameter:	0.50 mm
Inner Box(cm):	22x18x10.7
Unit:	10 Pcs
Export Carton:	10 Boxes
Export Carton Size(cm):	38x24x56



Skin Staple Remover

Description

The Skin Staple Remover is a sterile, ready-to-use medical tool designed for the efficient removal of surgical skin staples. Made from high-quality stainless steel, this device ensures durability and reliability in clinical settings. It is particularly recommended for use with the Universal Sutures Skin Stapler but is also compatible with most leading brands of skin staplers.

The Skin Staple Remover is suitable for use in various medical and surgical procedures where skin staples are employed for wound closure. Its design and material make it ideal for both initial and repeat use, provided proper sterilization protocols are followed.

Features

- Constructed from 100% pure stainless steel, providing strength and resistance to corrosion.
- Each Skin Staple Remover is sterile and ready for immediate use, ensuring patient safety and hygiene.
- The device is designed for ease of use, featuring an ergonomic shape that fits comfortably in the hand, reducing fatigue during prolonged use.
- The design allows for a clear view of the operating site, facilitating precise staple removal and minimizing the risk of tissue damage.
- The construction prevents the pins from breaking during use, ensuring consistent performance. It is compatible with various brands of skin staplers, making it a versatile tool in different clinical environments.




Absorbable Sutures

Description

Absorbable suture materials hold the wound edges in place during the healing phase. Absorbable materials support the natural healing process until form and function are restored. Such materials are subsequently metabolised by the organism. Absorbable suture materials are broken down either by endogenous proteolytic enzymes or by hydrolysis.

Features

- Knot security.
- High tensile strength.
- Problem-free disintegration and elimination.
- For tissue adaptations where the need for mechanical support is time-restricted.
- Optimum biocompatibility due to the chemical properties and refinements.



Absorbable Sutures

PARAMONO | Monofilament Poliglecaprone 25 U.S.P.
PARACRYL | Braided and Coated Polyglactin 910 U.S.P.
PARAQUICK | Braided & Coated Polyglactin 910- Fast Absorbable U.S.P.

PARASORB | Braided & Coated Polyglycolic Acid U.S.P.
PARADIOX | Monofilament Polydioxanone U.S.P.




Non-Absorbable Sutures

Description

Non-absorbable suture materials remain virtually unchanged in body tissues. Non-absorbable suture materials guarantee lasting support and best biotolerance, which is especially essential for long-term implants.

Features

- Long-term durability with simultaneous optimal biocompatibility.
- Optimal tissue compatibility in the body due to the careful selection of materials and refining process.
- Good tissue passage.



Non - Absorbable Sutures

PARALENE | Monofilament Polypropylene U.S.P.

PARALON | Braided & Coated Polyglycolic Acid U.S.P.

PARASILK | Black Braided Silk U.S.P.

PARABOND | Monofilament Polyamide (Nylon) U.S.P.



2 Way Standard Foley Catheter

Description

The Foley Catheter is a high-quality, medical-grade catheter designed for long-term urinary catheterization. Constructed from soft, bio compatible silicone, it provides exceptional comfort and durability for patients requiring continuous urinary drainage.

Features

- Smooth tapered tip facilities for easy introduction thereby reducing urethral trauma.
- Symmetrical shaped balloon which expands equally in all directions.
- Accurately formed drainage lumen for efficient flow of urine.
- Larger circular drainage lumen for efficient flow of urine.
- Smooth flexible siliconised catheter shaft to provide patient comfort.
- Integrally bonded Reinforced junction between the inflation funnel and the catheter shaft to ensure reliability thus preventing detachment during use.
- Plastic valves or Rubber valves fitted as per the convenience of the customers.

Details:


Product Name:	2 Way Foley Catheter
Unit Pack:	Box
Inner Box size (cm):	50x7.3x4.3
Inner Box:	10 Pcs
Export Carton size (cm):	44x44x51.5

Fitted with Hard Plastic Valve:

Size range:	12-30 FR/CH
-------------	-------------

Balloon Capacity range:

5-15ml/cc (12-22 FR)
15-30ml/cc (12-30 FR)
30ml/cc (12-30 FR)



2 Way Female Foley Catheter

Description

Similar in design to the 2 way standard foley catheters the latex female foleys are of a shorter length (Approximately 260mm) It is preferred for comfortable and convenient catheterisation of female patients.



Details:

Product Name:	2 Way Female Foley Catheter
Unit Pack:	Box
Inner Box size (cm):	50x7.3x4.3
Inner Box :	10 Pcs
Export Carton size (cm):	44x44x51.5

Fitted with Hard Plastic Valve:

Size range:	12 - 22 FR/CH
-------------	---------------

Balloon Capacity range:

5 - 15 ml/cc

3 Way Foley Catheter

Description

3 way catheters have been designed to irrigate the bladder continuously without disconnecting the catheter from a closed drainage system
3 way foley catheters are manufactured with highly advanced techniques to ensure

Features

- Maximum Patient Comfort
- Minimising infection thus increasing the confidence level of the customers.
- Saving of Valuable Nursing Time.
- Opposed Eyes of irrigation eye and drainage eye for efficient irrigation.
- Larger circular irrigation lumen for effective flow of irrigation fluids.
- 3 way foley catheters are adequately stiffer in order to withstand traction.

Details:

Product Name:	3 Way Foley Catheter
Unit Pack:	Box
Inner Box size (cm):	50x7.3x4.3
Inner Box:	10 Pcs
Export Carton size (cm):	44x44x51.5

Fitted with Hard Plastic Valve:

Size range: 14 - 26 FR/CH

Balloon Capacity range:

30 - 50 ml/cc



100% Silicone Foley Catheters - Latex Free

Description

100% Silicone Foley Catheters are made of highly biocompatible Silicone Elastomer, ideal for creating a safe, dependable Foley catheter for very long term implantation.

Features

- The non-tacky surface of Silicone makes it highly resistant to Urinary salt incrustation.
- Due to the inertness of Silicone, it does not leach chemicals, does not irritate Urethral Mucosa and also does not lead to Urethral stricture formation.
- The smooth and longer eyes allow non-traumatic insertion and maximum drainage.
- The ribbed funnel ensures an easy, leak-proof connection while the clear connector and shaft allow visual examination of drainage.
- The drainage lumen is comparatively bigger and thus allows better drainage of Urine and blood clots.

Our Product Ranges - 100% Silicone Foley Catheters

Type	Shaft Gauge	Balloon Capacity (cc/ml)
2way Paediatric Special	6 fr	1.5 cc/ml
2way Paediatric	8 to 10 fr	3 cc/ml
2way Standard	12 to 24 fr	5 - 15 cc/ml
2way Standard	16 to 26 fr	30 cc/ml
3way Standard	16 to 24 fr	30 cc/ml



Microtome Blade


Description

Precision, high-quality workmanship, superior service life – the blades from Paramount offer the optimum conditions for delivering the best work in the medical field.

Paramount | PN35


The durable blade for enduring cuts in hard tissue

For thin sectioning of paraffin-embedded tissue blocks. The blade is characterized by high durability.



Details:

Product Name:	Microtome Blade
Inner Box(cm):	25X12.5X5.8
Unit:	10Pcs
Export Carton:	10 Boxes
Export Carton Size(cm):	27X27X31



Paramount | PA35

The sharp blade for precise cuts in hard and soft tissue

For very fine cuts of paraffin-embedded tissue blocks due to high sharpness. Long service life even when cutting through hard tissue.



Paramount | PS35

The universal blade for stable cutting quality

For thin sectioning of paraffin-embedded tissue blocks. The blade features superior cutting performance and durability.



Blood Lancet

Description

Supplied sterile and individually packed. Lancet comes in a pack of 200 pieces per box. the tip is specially designed for easy pricking, which minimizes patient trauma.



Details:

Product Name:	Blood Lancet
Product Code:	PSSBL
Size of Inner Box (cm.):	12.5 x 7 x 20
Unit/Box:	200
Export carton:	100 Boxes x 200 Units
Export Carton size (cm.):	44.5 x 38.5 x 14
Per 20 feet container:	1800 Cartons



Blood Lancet - Twisted

Description

Supplied sterile and individually packed. Lancet comes in a pack of 100 pieces per box. the tip is specially designed for easy pricking, which minimizes patient trauma.



Details:

Product Name:	Blood Lancet - Twisted
Product Code:	PSSBLT
Size of Inner Box (cm.):	6.5 x 4.2 x 7.0
Unit/Box:	100
Export carton:	64 Boxes x 100 Units
Export Carton size (cm.):	30 x 27.5 x 17
Per 20 feet container:	1800 Cartons

Safety Lancet

Features

- Single use only
- Sterilized by Gamma Radiation 2.5 M Rads
- With 30G Needle



Flyfix

Features

- Discard after single use.
- Convenient.
- Extra strip for strong adhesion
- Saves nursing time.
- Infusion plaster with wound pad.



Details:

Product Name:	I.V Fixator
Product Code:	PAFF-24
Size of Box (cm.):	12x11x10
Unit/Box:	50
Export carton:	30 Boxes x 50 Units
Export carton size (cm):	51 x 39 x 25

Disposable Adult Diapers


For Urinary Incontinence Management

Features:

- Dry-Fast Core, Super absorbent.
- Moisture repellent.
- Security barriers.
- Waist Elastication and curved leg elastic.
- Superfit Band Wetness Indicator.
- Three absorbency levels and three sizes.
- Compact and convenient packaging.
- Also available in options of double cotton core & ADL.

Benefits:

- Faster absorption, drier surface.
- Better containment of faces.
- Reduces urine leakage.
- Improves comfort.
- Comfortable fit.
- Easier to fix and adjust tapes.
- A Guide for carers.
- Meets individual requirements with cost effective care.
- Easy to handle less storage space required.



Details:

Product Name:	Adult Diapers (Medium)
Product Code:	AKMP 04
Size of Poly Bag (cm.):	20 x 23.0 x 17.5
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	35.5x 46 x 64
Per 20 feet container:	240 Cartons

Details:

Product Name:	Adult Diapers (Large)
Product Code:	AKLP 06L
Size of Poly Bag (cm.):	22 x 24 x 17.5
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	48.5 x 35.5 x 68
Per 20 feet container:	220 Cartons

Details:

Product Name:	Adult Diapers (Extra Large)
Product Code:	AKXLP 07
Size of Poly Bag (cm.):	23.5 x 24 x 17.5
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	49 x 35.5 x 74
Per 20 feet container:	200 Cartons

**INDIA'S FIRST
ISI APPROVED
ADULT DIAPERS**

IS : 17508:2020



CM/L 8600113014




Disposable Adult Diapers

Premium



ADL (Acquisition Distribution Layer)



Details:

Product Name:	Pullups Medium
Product Code:	APM10
Size of Poly Bag (cm.):	15X18.5X18.5
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	51X47X36.5
Per 20 feet container:	290 Cartons

Economy



Details:

Product Name:	Pullups Large
Product Code:	APL11
Size of Poly Bag (cm.):	16X19X21
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	53X46.5X41
Per 20 feet container:	250 Cartons

Details:

Product Name:	Pullups X-Large
Product Code:	APXL12
Size of Poly Bag (cm.):	20X16X21.5
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size (cm.):	54X48.5X42
Per 20 feet container:	230 Cartons



Under Pads

AroKleen

Description

- **Anti-Leakage Design:**
The product has extra width and length designed for better absorption.
- **Dry And Comfortable:**
It has large absorbing padding using the reticular embossed design to introduce the water diversion towards every side keep the surface comfortable.
- **Skin Friendly:**
Using the super-nonwoven as the top layer, it feels comfortable and prevents skin rashes.
- **International Quality:**
Strict controls of product material and production assuring safety and health.

Details:

Product Name:	AroKleen Baby Mat
Product Code:	ABM60x60
Size of Poly Bag (cm.):	29 x 20 x 10
Unit/Pkt:	10
Export carton:	15 Pkt x 10 Units
Export Carton size (cm.):	37 x 28 x 56
Per 20 feet container:	440 Cartons


Details:

Product Name:	Under Sheet
Product Code:	US90x60
Size of Poly Bag (cm.):	29 x 20 x 10
Unit/Pkt:	10
Export carton:	12 Pkt x 10 Units
Export Carton size:	41 x 27.5 x 55
Per 20 feet container:	410 Cartons

Details:

Product Name:	Puppy Pad
Product Code:	PP30P
Size of Poly Bag (cm.):	26.5 x 18 x 27
Unit/Pkt:	30
Export carton:	4 Pkt x 30 Units
Export Carton size (cm.):	60 x 36 x 29
Per 20 feet container:	290 Cartons

5 LAYERED LEAK PROOF TECHNOLOGY!



The Advantages of Using PupPee Pads

Our pads provide worry-free training for your puppy & peace of mind for you. They also can be used as an alternative to outdoor "relief" for your adult dog.

- Locks in wetness with quilted top layer
- Ensures maximum absorption power with 2 moisture-locking layers
- Protects floor & carpets with leak-proof backing
- Provides easy training with built-in attractant
- Protect your floors
- Easy cleanup & disposal
- Quilted for extra absorbency
- Less accidents means less disciplining
- Dog attractant reinforces training
- Durable, tear resistant



Disposable Prep Razor

Paraprep

Features

- Safe contoured handle with sure grip ridges provide precise fingertips control.
- Easy to handle lies in front of blades thus requiring only light gentle pull instead of pressure.
- Coated stainless steel blades shaves close without nick or irritation.
- Open design provides optimal view of shaving area.
- Compact and convenient packaging.
- Cost effective, no need for more than one razor per procedure.

Details:

Product Name:	Prep Razor
Product Code:	PPDR20210320
Size of Poly Bag (Cm.):	9.5x5.7
Inner Box:	20x12x6.5
Unit/Box:	50
Export carton:	20 Box x 50 Units
Export Carton size:	41.5X27.5X35
Per 20 feet container:	650 Cartons



with Two Side Blades



For All Body
Areas



Smooth Hair
Removal



100% Shower
Proof



Global Standard
Manufacturing



Unique Skin
Protection





Stainless Steel
Blade

Details:

Product Name:	2 Side Prep Razor
Product Code:	PPRD
Size of Poly Bag (Cm.):	8.5x15.5
Inner Box:	28.5x19X26
Unit/Box:	100
Export carton:	5 Box x 100 Units
Export Carton size:	59 X 49 X 28
Per 20 feet container:	650 Cartons



Face Mask - 3 Ply/ Type IIR



3 PLY Mask - Non Surgical



Smooth Fiber
Material



SMS Fabric

3 PLY Mask - Surgical



Meltblown
Fabric



UV Cleaned

FEATURES



Adjustable Nose
Clip



Ultrasonic
Sealing



One Size
Fits Most



Seamless
Crimping



Comfortable
& Breathable



Quality
Material



3 Layer
Protection



Anti Droplets

Details:

Product Name:	3 PLY/ Type IIR Mask - Non Surgical
Product Code	P03N
Size of Inner Box (cm.):	23 x 11 x 19.5
Unit/Box:	100
Export Carton size (cm.):	57 x 48 x 41
Per 20 feet container:	230 Carton

Details:

Product Name:	3 PLY/ Type IIR Mask - Surgical
Product Code	P03M
Size of Inner Box (cm.):	23 x 10.5 x 10
Unit/Box:	50
Export Carton size (cm.):	55 x 48 x 42
Per 20 feet container:	230 Carton

Face Mask - N95/FFP2



without Respirator

with earloop



UV Cleaned



without Respirator

with headloop



Filters Particulates



BFE≥95%



with Respirator

with earloop



PFE≥95%

FEATURES



Adjustable Noise Clip




Soft Elastic Ear Loops



High Quality Breathable
Material and Valve



Seamless Flow of
Exhaled Air



Details:

Product Name:	N95/FFP2 without Respirator
Product Code:	P95W
Size of Inner Box (cm.):	30.5 x 13.5 x 17.5
Unit/Box:	50
Export Carton size (cm.):	63.5 x 56 x 55
Per 20 feet container:	140 Carton

Details:

Product Name:	N95/FFP2 with Respirator
Product Code:	P95R
Size of Inner Box (cm.):	24.5 x 13.5 x 17.5
Unit/Box:	25
Export Carton size (cm.):	51 x 42 x 54.5
Per 20 feet container:	220 Carton

Bed Bath Wipes

AroKleen

Description

AroKleen Bed Bath Wipes are Paraben Free, Alcohol Free, SLS Free, Microbicide. The fast acting formula, with a two minute kill time loosens and suspends blood and bodily fluids, non-organic soil and other cellular debris for easy cleaning.

CLEANS, DISINFECTS & DEODORIZES

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



Skin Friendly
pH



Aloe Vera



Enriched
With Glycerin

USAGES



Cleansing



Refreshing



Moisturizing



Hospital Care

Cleansing & Refreshing Wipes

SPECIALLY FOR CLEANSING

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



No Added
Alcohol



Skin Friendly
pH



Enriched With
Vitamin C



Details:

Product Name:	Bed Bath Wipes
Product Code:	BP2202
Size of Poly Bag (Cm.):	10X8.5X24
Unit/Pkt:	10
Export carton:	48 Pkt x 10 Units
Export Carton size(cm):	24x22x44
Per 20 feet container:	1100 Cartons



**Enriched With
2% Chlorhexidine**

Details:

Product Name:	Cleansing Wipes
Product Code:	BN2202
Size of Poly Bag (Cm.):	10X9.2X23.5
Unit/Pkt:	30
Export carton:	48 Pkt x 30 Units
Export Carton size(cm):	36 x 32 x 22
Per 20 feet container:	1000 Cartons

Germ Protection Wipes

AroKleen

Description

AroKleen Germ Protection Wipes are more effective and cost less than sprays and towels, making them the right wipe to use in any location. Every area where staff, visitors, guests, customers and the general public convene, there are surfaces that need to be regularly cleaned and freed of dirt, grime, harmful germs and bacteria.


PROTECTION FROM GERMS

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



USAGES



Details:

Product Name:	Germ Protection Wipes
Product Code:	PGP01
Size of Poly Bag (Cm.):	10X8.5X24
Unit/Pkt:	48
Export carton:	36 Pkt x 48 Units
Export Carton size(cm):	36 x 33 x 23
Per 20 feet container:	920 Cartons



Baby Wipes

AroKleen

Description

AroKleen Baby Wipes are large, durable, disposable clothes that are alcohol and latex free. Our Aloe Vera Extract, Vitamin E lotion-based formula that gently cleans and keeps babies' skin smooth and moisturized.

SPECIALLY FOR BABIES

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



Details:

Product Name:	Baby Wipes
Product Code:	BB23009
Size of Poly Bag (Cm.):	10X10.3X24.5
Unit/Pkt:	72
Export carton:	24 Pkt x 72Units
Export Carton size (cm):	37 x 32 x 25
Per 20 feet container:	850 Cartons



Baby Wipes with 99.9 % Water

SPECIALLY FOR KIDS

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



No Added
Parabens



Aloe Vera



Silicon Free



pH Balanced



No Synthetic
Preservatives

Details:

Product Name:	Baby Wipes
Product Code:	BWW01
Size of Poly Bag (Cm.):	10X10X24
Unit/Pkt:	72
Export carton:	24 Pkt x 72 Units
Export Carton size:	36 x 32 x 23.5
Per 20 feet container:	850 Cartons



Pet Wipes

SPECIALLY FOR CATS & DOGS

Quick | Easy | Effective | Kills Virus

UNIQUE FEATURES



Powerful Odour



pH Balanced



Olive Oil



Enriched With
Vitamin E



Cruelty Free



Vet & Groomer
Recommended

Details:

Product Name:	Pet Wipes
Product Code:	PGP01
Size of Poly Bag (Cm.):	10X10X24
Unit/Pkt:	48
Export carton:	36 Pkt x 48 Units
Export Carton size:	38 x 32 x 22
Per 20 feet container:	950 Cartons



Lung Exerciser

Spirometer

Features:

- Three color balls for easier identification of inhale capacity.
- Improve maximal exercise capacity.
- Made from transparent material.
- Compact, convenient and safe.

Details:

Product Name:	Lung Exerciser
Product Code:	PLE01
Size of Box (Cm.):	14.2x7x15.7
Unit/Box:	01
Export carton:	12 Bags x 01 Units



Hot Water Bottle



Joint Pain



Sports Pain



Arthritis Pain



Back Aches



Stomach Aches



Muscular Cramps



Rheumatic Pain



Muscle Pull



Regular Hot Water
Bottle

Product Code: PHWB 05



BS Standard Hot
Water Bottle

Product Code: HWBP 06

Infrared Forehead Thermometer

Features

- Non-contact Forehead
- LCD Backlight
- C/F Unit Switchable

Details:

Product Name:	Infrared Forehead Thermometer
Product Code:	PITM-10
Size of Box (Cm.):	10.5X4.65X16.5
Unit/Box:	01
Export carton:	12 Bags x 01 Units



16 Memories



1s Fast Reading



High Temperature Alarm



Pulse Oximeter

Features

- 2pcs standard AAA1.5V alkaline battery supporting 20 hours continuous operation.
- Applicable to Hospital, Clinic, Family, Oxygen Bar, Sporting center etc.
- True color high resolution OLED/TFT LCD/LED.
- Light in weight, simple in operation, easy to use & carry.
- Adjust the parameters in friendly menu.
- Soft anti-ambient light interference.
- Low power consumption.
- Optional Pouch.

Details:

Product Name:	Pulse Oximeter
Product Code:	PPOM-09
Size of Box (Cm.):	8.5x5x6
Unit/Box:	01
Export carton:	12 Bags x 01 Units



Air Mattress

Bedsore Prevention Cure Air Mattress

Features

- Anti decubitus alternating pressure air mattress Provide relief for bedsore.
- Low noise, low vibration & low power Consumption.
- Recommended for bedridden patients.
- Medical grade PVC.
- Long working life.
- Large air output.

Details:

Product Name:	Air Mattress
Product Code:	PACM-24
Size of Box(cm.):	26x16x20
Unit/Box:	01
Export carton:	06 Bags x 01 Units
Export carton Size:	50 x 28 x 42.5



16 Memories



High Temperature Alarm



1s Fast Reading





1, L.S.C., Okhla Industrial Area, Phase-II, New Delhi-110 020, India

Tel +91-11-46436600, 46436601, 46436666, 26389812, 13

Fax +91-11-26389815, 41616555

sales@paramountblades.com

export2@paramountblades.com

www.paramountsurgimed.com

(Disclaimer- All measurements in weight & size are indicative for reference only.)

PSL/CTLG/01, Issue No. 01, Issue Date : 26/11/2020, Rev. No. 05, Rev. Date : 24/03/2025